

Angle-stable distal radius plates

for screws Ø 3.5 mm; "Quadro"; titanium



Distal radius plate 3.5 "Quadro"
Item no. 5.876.43

Item no. Item description

Plate thickness 1.7 mm; conical thread in head

5.876.43	Distal radius plate 3.5 "Quadro"; 4 head holes; 3 shaft holes; 51 mm x 10 mm; right
5.876.44	Distal radius plate 3.5 "Quadro"; 4 head holes; 4 shaft holes; 59 mm x 10 mm; right
5.877.43	Distal radius plate 3.5 "Quadro"; 4 head holes; 3 shaft holes; 51 mm x 10 mm; left
5.877.44	Distal radius plate 3.5 "Quadro"; 4 head holes; 4 shaft holes; 59 mm x 10 mm; left

Plate thickness 2.0 mm; conical thread in head and shaft

5.876.44/2	Distal radius plate 3.5 "Quadro"; 4 head holes; 4 shaft holes; 59 mm x 10 mm; right
5.876.46/2	Distal radius plate 3.5 "Quadro"; 4 head holes; 6 shaft holes; 76 mm x 10 mm; right
5.876.48/2	Distal radius plate 3.5 "Quadro"; 4 head holes; 8 shaft holes; 93 mm x 10 mm; right
5.877.46/2	Distal radius plate 3.5 "Quadro"; 4 head holes; 6 shaft holes; 76 mm x 10 mm; left
5.877.48/2	Distal radius plate 3.5 "Quadro"; 4 head holes; 8 shaft holes; 93 mm x 10 mm; left

on request:

5.876.43/2	Distal radius plate 3.5 "Quadro"; 4 head holes; 3 shaft holes; 51 mm x 10 mm; right
5.877.44/2	Distal radius plate 3.5 "Quadro"; 4 head holes; 4 shaft holes; 59 mm x 10 mm; left

Plate thickness 3.0 mm; conical thread in head and shaft

on request:

5.876.48/3	Distal radius plate 3.5 "Quadro"; 4 head holes; 8 shaft holes; 93 mm x 10 mm; right
5.877.48/3	Distal radius plate 3.5 "Quadro"; 4 head holes; 8 shaft holes; 93 mm x 10 mm; left

Information:

Combi holes in the shaft for optional fitting of angle-stable screws or standard screws.



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498 - 560
Fax +49 36738 498 - 89 560

Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com