

Variable Angle-stable Forefoot Step-AD Plates

Small Fragment; for Screws \varnothing 2.7 mm; Titanium



Item no. 5.2913.22/2



Item no. 5.2913.22/3



Item no. 5.2913.22/4



Item no. 5.2913.22/5

Item no.	Item description
5.2913.22/2	step-AD plate with conical thread \varnothing 2.7 mm; 5 holes; 22 mm x 16 mm; offset 2 mm; variable angle-stable; titanium; gold
5.2913.22/3	step-AD plate with conical thread \varnothing 2.7 mm; 5 holes; 22 mm x 16 mm; offset 3 mm; variable angle-stable; titanium; pink
5.2913.22/4	step-AD plate with conical thread \varnothing 2.7 mm; 5 holes; 22 mm x 16 mm; offset 4 mm; variable angle-stable; titanium; light blue
5.2913.22/5	step-AD plate with conical thread \varnothing 2.7 mm; 5 holes; 22 mm x 16 mm; offset 5 mm; variable angle-stable; titanium; dark purple

Information:

Possible combination of these plates with the set variable angle-stable forefoot plates 2.7, set no. 19.523.00, titanium.
Central non angle-stable hole for compression.