

OP-Set Dens-Axis - screwing



■ Fracture compressing screw, Titanium

Diameter ; length	Code-No.	Quantity
Diameter 3.0/4.0 mm; L 32 mm	3.123.32	3
Diameter 3.0/4.0 mm; L 34 mm	3.123.34	3
Diameter 3.0/4.0 mm; L 36 mm	3.123.36	3
Diameter 3.0/4.0 mm; L 38 mm	3.123.38	3
Diameter 3.0/4.0 mm; L 40 mm	3.123.40	3
Diameter 3.0/4.0 mm; L 42 mm	3.123.42	3
Diameter 3.0/4.0 mm; L 44 mm	3.123.44	3



Bestell-Nr.
Komplettsset
Titan

19.380.10

■ Instruments

Item designation	Code No.	Quantity
Cannulated pilot drill with handle and cutter, D 2.1 x 352 mm	10.266.01	1
Cannulated hex screwdriver with handle, wrench width 2.0 mm, length 350 mm	10.266.02	1
Double tissue protecting drill sleeve D 7.0 x 140 mm	10.266.03	1
Sleeve for 10.266.03	10.266.034	2
Gauge, cannulated	10.266.04	1
Screw holding forceps	2.954.01	1
Tray for implants and instruments for dens-screwing	19.381.10	1



■ Guidewire

Item designation	Code No.	Quantity
Guidewire, thread 10 mm diameter 1.2 mm, L 400 mm	2.984.12/400	5

OP-Set Dens-Axis - screwing

