


KI-CVA-Platten

für Schrauben Ø 4.5/5.2 mm; Titan



Art.-Nr. 5.555.1624



Art.-Nr. 5.555.1640



Art.-Nr. 5.555.1652


Artikel-Nr.	Bezeichnung	EAN
5.555.1624	KI-CVA Platte 4.5/5.2; 4 Loch; 16 x 24 mm	
5.555.1626	KI-CVA Platte 4.5/5.2; 4 Loch; 16 x 26 mm	
5.555.1628	KI-CVA Platte 4.5/5.2; 4 Loch; 16 x 28 mm	
5.555.1630	KI-CVA Platte 4.5/5.2; 6 Loch; 16 x 30 mm	
5.555.1632	KI-CVA Platte 4.5/5.2; 6 Loch; 16 x 32 mm	
5.555.1634	KI-CVA Platte 4.5/5.2; 6 Loch; 16 x 34 mm	
5.555.1636	KI-CVA Platte 4.5/5.2; 6 Loch; 16 x 36 mm	
5.555.1638	KI-CVA Platte 4.5/5.2; 6 Loch; 16 x 38 mm	
5.555.1640	KI-CVA Platte 4.5/5.2; 6 Loch; 16 x 40 mm	
5.555.1642	KI-CVA Platte 4.5/5.2; 6 Loch; 16 x 42 mm	
5.555.1644	KI-CVA Platte 4.5/5.2; 6 Loch; 16 x 44 mm	
5.555.1646	KI-CVA Platte 4.5/5.2; 8 Loch; 16 x 46 mm	
5.555.1648	KI-CVA Platte 4.5/5.2; 8 Loch; 16 x 48 mm	
5.555.1650	KI-CVA Platte 4.5/5.2; 8 Loch; 16 x 50 mm	
5.555.1652	KI-CVA Platte 4.5/5.2; 8 Loch; 16 x 52 mm	
5.555.1654	KI-CVA Platte 4.5/5.2; 8 Loch; 16 x 54 mm	
5.555.1657	KI-CVA Platte 4.5/5.2; 8 Loch; 16 x 57 mm	
5.555.1660	KI-CVA Platte 4.5/5.2; 8 Loch; 16 x 60 mm	
5.555.1663	KI-CVA Platte 4.5/5.2; 10 Loch; 16 x 63 mm	
5.555.1666	KI-CVA Platte 4.5/5.2; 10 Loch; 16 x 66 mm	
5.555.1669	KI-CVA Platte 4.5/5.2; 10 Loch; 16 x 69 mm	
5.555.1672	KI-CVA Platte 4.5/5.2; 10 Loch; 16 x 72 mm	
5.555.1675	KI-CVA Platte 4.5/5.2; 10 Loch; 16 x 75 mm	

KI-CVA-Platten

für Schrauben Ø 4.5/5.2 mm; Titan









Art.-Nr. 5.555.1624



Art.-Nr. 5.555.1640



Art.-Nr. 5.555.1652

Artikel-Nr.	Bezeichnung	EAN
5.555.1678	KI-CVA Platte 4.5/5.2; 12 Loch; 16 x 78 mm	 4 047075 164042
5.555.1681	KI-CVA Platte 4.5/5.2; 12 Loch; 16 x 81 mm	 4 047075 164059
5.555.1684	KI-CVA Platte 4.5/5.2; 12 Loch; 16 x 84 mm	 4 047075 164066
5.555.1687	KI-CVA Platte 4.5/5.2; 12 Loch; 16 x 87 mm	 4 047075 164073
5.555.1690	KI-CVA Platte 4.5/5.2; 12 Loch; 16 x 90 mm	 4 047075 164080
5.555.1693	KI-CVA Platte 4.5/5.2; 14 Loch; 16 x 93 mm	 4 047075 164097