

Distale Radiusplatte 2.4 vws Watershed

variabel winkelstabil; für Schrauben Ø 2.4 mm; Titan



Art.-Nr. 5.08373.03







Art.-Nr. 5.08373.04



Art.-Nr. 5.08473.03



Art.-Nr. 5.08473.04

Artikel-Nr.	Bezeichnung	EAN
5.08373.03	Dist.Radiuspl. Watershed 2.4; 7/3 Loch; 23x50; schmal; vws; links	 4 047075 189489
5.08373.04	Dist.Radiuspl. Watershed 2.4; 7/4 Loch; 23x59; schmal; vws; links	 4 047075 189496
5.08473.03	Dist.Radiuspl. Watershed 2.4; 7/3 Loch; 23x50; schmal; vws; rechts	 4 047075 189526
5.08473.04	Dist.Radiuspl. Watershed 2.4; 7/4 Loch; 23x59; schmal; vws; rechts	 4 047075 189533

Information

Sie können die Platte mit dem Instrumentarium folgender Sets verwenden: 19.117.00 Set Distale Radiusplatte 2.4 vws und 19.024.001 Set Instrumente 2.4

