

# Distale Radiusplatte 3.5 vws

variabel winkelstabil; für Schrauben Ø 3.5 mm; Titan



Art.-Nr. 5.8763.34

Art.-Nr. 5.8823.34

Artikel-Nr.	Bezeichnung	EAN
-------------	-------------	-----

### Distale Radiusplatte 3.5 Classic vws

5.8763.33	Dist. Radiusplatte Classic 3.5; 3/3 Loch; 51x10mm; R; Stärke 3,0mm Kopf/2,0mm Schaft	 4047075111428
5.8763.34	Dist. Radiusplatte Classic 3.5; 3/4 Loch; 59x10mm; R; Stärke 3,0mm Kopf/2,0mm Schaft	 4047075111435
5.8763.36	Dist. Radiusplatte Classic 3.5; 3/6 Loch; 76x10mm; R; Stärke 3,0mm Kopf/2,0mm Schaft	 4047075121106
5.8763.38	Dist. Radiusplatte Classic 3.5; 3/8 Loch; 93x10mm; R; Stärke 3,0mm Kopf/3,0mm Schaft	 4047075121113
5.8763.310	Dist. Radiusplatte Classic 3.5; 3/10 Loch; 110x10mm; R; Stärke 3,0mm Kopf/3,0mm Schaft	 4047075121083
5.8763.313	Dist. Radiusplatte Classic 3.5; 3/13 Loch; 136x10mm; R; Stärke 3,0mm Kopf/3,0mm Schaft	 4047075121090
5.8773.33	Dist. Radiusplatte Classic 3.5; 3/3 Loch; 51x10mm; L; Stärke 3,0mm Kopf/2,0mm Schaft	 4047075111442
5.8773.34	Dist. Radiusplatte Classic 3.5; 3/4 Loch; 59x10mm; L; Stärke 3,0mm Kopf/2,0mm Schaft	 4047075111459
5.8773.36	Dist. Radiusplatte Classic 3.5; 3/6 Loch; 76x10mm; L; Stärke 3,0mm Kopf/2,0mm Schaft	 4047075121144
5.8773.38	Dist. Radiusplatte Classic 3.5; 3/8 Loch; 93x10mm; L; Stärke 3,0mm Kopf/3,0mm Schaft	 4047075121151
5.8773.310	Dist. Radiusplatte Classic 3.5; 3/10 Loch; 110x10mm; L; Stärke 3,0mm Kopf/3,0mm Schaft	 4047075121120
5.8773.313	Dist. Radiusplatte Classic 3.5; 3/13 Loch; 136x10mm; L; Stärke 3,0mm Kopf/3,0mm Schaft	 4047075121137

### Distale Radiusplatte 3.5 Minor vws

5.8823.33	Dist. Radiusplatte Minor 3.5; 3/3 Loch; 51x10mm; R; Stärke 3,0mm Kopf/2,0mm Schaft	 4047075116492
5.8823.34	Dist. Radiusplatte Minor 3.5; 3/4 Loch; 59x10mm; R; Stärke 3,0mm Kopf/2,0mm Schaft	 4047075116508
5.8823.36	Dist. Radiusplatte Minor 3.5; 3/6 Loch; 76x10mm; R; Stärke 3,0mm Kopf/2,0mm Schaft	 4047075120338
5.8823.38	Dist. Radiusplatte Minor 3.5; 3/8 Loch; 93x10mm; R; Stärke 3,0mm Kopf/3,0mm Schaft	 4047075120345
5.8823.310	Dist. Radiusplatte Minor 3.5; 3/10 Loch; 110x10mm; R; Stärke 3,0mm Kopf/3,0mm Schaft	 4047075120314
5.8823.313	Dist. Radiusplatte Minor 3.5; 3/13 Loch; 136x10mm; R; Stärke 3,0mm Kopf/3,0mm Schaft	 4047075120321
5.8833.33	Dist. Radiusplatte Minor 3.5; 3/3 Loch; 51x10mm; L; Stärke 3,0mm Kopf/2,0mm Schaft	 4047075116515
5.8833.34	Dist. Radiusplatte Minor 3.5; 3/4 Loch; 59x10mm; L; Stärke 3,0mm Kopf/2,0mm Schaft	 4047075116522
5.8833.36	Dist. Radiusplatte Minor 3.5; 3/6 Loch; 76x10mm; L; Stärke 3,0mm Kopf/2,0mm Schaft	 4047075120376
5.8833.38	Dist. Radiusplatte Minor 3.5; 3/8 Loch; 93x10mm; L; Stärke 3,0mm Kopf/3,0mm Schaft	 4047075120383
5.8833.310	Dist. Radiusplatte Minor 3.5; 3/10 Loch; 110x10mm; L; Stärke 3,0mm Kopf/3,0mm Schaft	 4047075120352
5.8833.313	Dist. Radiusplatte Minor 3.5; 3/13 Loch; 136x10mm; L; Stärke 3,0mm Kopf/3,0mm Schaft	 4047075120369