

Set CPL plate 4.5 as IH


Cable-Plate-System Labitzke; angle-stable; stainless steel; Set no. 19.516.00



Item no. 4.2341.10





Item no. 4.2461.10



Set no. 19.516.00  Internal hexagon

Item no.	Item description	Quant.
Implants – plates		
CPL plate with conical thread 4.5; stainless steel		
4.2341.08	6 notches; 8 holes; 124 x 12 mm	1
4.2341.10	8 notches; 10 holes; 154 x 12 mm	1
4.2341.12	10 notches; 12 holes; 184 x 12 mm	1
4.2341.14	12 notches; 14 holes; 214 x 12 mm	1
4.2461.10	8 notches; 10 holes; 174 x 16 mm	1
4.2461.12	10 notches; 12 holes; 206 x 16 mm	1
4.2461.14	12 notches; 14 holes; 238 x 16 mm	1
4.2461.16	14 notches; 16 holes; 270 x 16 mm	1

Item no.	Item description	Quant.
Implants – screws		
Cortical screw Ø 4.5 mm; fully threaded; self tapping; stainless steel		
 Internal hexagon		
2.152.12	length 12 mm	3
2.152.14	length 14 mm	3
2.152.16	length 16 mm	3
2.152.22	length 22 mm	4
2.152.24	length 24 mm	4
2.152.26	length 26 mm	4
2.152.28	length 28 mm	4
2.152.30	length 30 mm	4
2.152.32	length 32 mm	4
2.152.34	length 34 mm	4
2.152.36	length 36 mm	4
2.152.38	length 38 mm	4
2.152.40	length 40 mm	4
2.152.42	length 42 mm	4


Item no.	Item description	Quant.
Implants – screws		
Cortical screw Ø 4.5 mm; conical head thread; fully threaded; self tapping; stainless steel		
 Internal hexagon		
2.545.18	length 18 mm	5
2.545.20	length 20 mm	5
2.545.22	length 22 mm	5
2.545.24	length 24 mm	5
2.545.26	length 26 mm	5
2.545.28	length 28 mm	5
2.545.30	length 30 mm	5
2.545.32	length 32 mm	5
2.545.34	length 34 mm	5
2.545.36	length 36 mm	5
2.545.38	length 38 mm	5
2.545.40	length 40 mm	5
2.545.42	length 42 mm	5
2.545.45	length 45 mm	5
2.545.50	length 50 mm	5
2.545.55	length 55 mm	5

Set CPL plate 4.5 as IH

Cable-Plate-System Labitzke; angle-stable; stainless steel; Set no. 19.516.00

Item no.	Item description	Quant.
Instruments		
2.904.09	Drill bit qc 3.2; Spiral 60; 2-flute; length 195 mm	1
2.977.02	drill guide for angle stable screwing, length 70 mm, big fragment	2
10.500.21	drill guide for neutral and compression position, Ø 3.2 mm	1
10.502.32	universal drill guide, Ø 3.2 mm	1
2.953.90	gauge with clasp for screws Ø 4.5 mm, Ø 6.0 mm, measuring range 90 mm	1
2.9406.35	hex screwdriver with handle for screws Ø 4.5 mm - Ø 7.0 mm, AF 3.5, length 215 mm	1
2.954.01	screw forceps, self holding	1

Item no.	Item description	Quant.
Container		
19.515.00	perforated autoclavable container with inset for instruments and CPL-plates	1



Set no. 19.516.00  Internal hexagon