



Set Calcaneus Plate 3.5/4.0 vas


variabel angle-stable; for screws 3.5/4.0 mm; titanium; Set no. 19.326.00




Item no. 5.894.00




Item no. 5.895.00



Item no. 5.8783.09



Item no. 5.8793.09



Set no. 19.326.00 Internal hexagon

Item no.	Item description	Quant.
Implants – plates		
Calcaneus plate 3.5/4.0 vas with extension; titanium		
5.894.00	calcaneus plate w. conic thread 3,5 a. 4,0; 10 holes; length 72 mm; vas, right	1
5.895.00	calcaneus plate w. conic thread 3,5 a. 4,0; 10 holes; length 72 mm; vas, left	1
Calcaneus plate 3.5/4.0 vas with bridges; titanium		
5.8783.09	calcaneus plate w. conic thread 3,5 a. 4,0; 9 holes; length 65 mm; vas, right	1
5.8793.09	calcaneus plate w. conic thread 3.5 a. 4.0, 9 holes, length 65 mm, vas, left	1


Item no.	Item description	Quant.
Implants – screws		
Cannulated screw Ø 4.0 mm; thread 14; self drilling, self tapping; back tapping; titanium		
3.440.40	length 40 mm	3
3.440.42	length 42 mm	3
3.440.44	length 44 mm	3
3.440.46	length 46 mm	3
3.440.48	length 48 mm	3
3.440.50	length 50 mm	3
3.440.55	length 55 mm	3

Implants – screws		
Cancellous screw Ø 4.0 mm; fully threaded; titanium		
3.145.22	length 22 mm	3
3.145.26	length 26 mm	3
3.145.30	length 30 mm	3
3.145.35	length 35 mm	3
3.145.40	length 40 mm	3
3.145.45	length 45 mm	3
3.145.50	length 50 mm	3

Cortical screw with conical head thread Ø 3.5 mm; fully threaded; self tapping; titanium		
3.133.22	length 22 mm	5
3.133.26	length 26 mm	5
3.133.30	length 30 mm	5
3.133.35	length 35 mm	5
3.133.40	length 40 mm	5
3.133.45	length 45 mm	5
3.133.50	length 50 mm	5

Set Calcaneus Plate 3.5/4.0 vas

variabel angle-stable; for screws 3.5/4.0 mm; titanium; Set no. 19.326.00



Set no. 19.326.00 Internal hexagon

Item no.	Item description	Quant.
Implants – screws		
	Cancellous screw with conical head thread Ø 4.0 mm; fully threaded; self tapping; titanium Internal hexagon	
3.147.22	length 22 mm	5
3.147.26	length 26 mm	5
3.147.30	length 30 mm	5
3.147.35	length 35 mm	5
3.147.40	length 40 mm	5
3.147.45	length 45 mm	5
3.147.50	length 50 mm	5

Wires		
6.031.16	Kirschner wire with trocar point and round end, Ø 1.6 x 150 mm, stainless steel	5

Container		
19.325.00	perforated autoclavable container with inset for instruments/implants for calcaneus-set	1

Item no.	Item description	Quant.
Instruments		
2.950.150	gauge for cannulated screws Ø 2.7 - 4.0 mm, L 150 mm	1
2.9436.40	hex screwdriver w. handle for Ø 2.7 - 4.0 mm; cannulated; I Ø 1.7 mm; wrench size 2.5 mm; L 215 mm	1
2.9406.25	hex screwdriver with handle for screws Ø 2.7 - 4.0 mm, AF 2.5, conical	1
2.953.60	gauge with clasp for screws Ø 3.5 mm, Ø 4.0 mm, measuring range 60 mm	1
2.977.08	drill guide w. handle f. variable angle stable screwing, small fragment, length 45 mm	1
2.977.01	drill guide for angle stable screwing, length 45 mm, small fragment	2
2.954.01	screw forceps, self holding	1
2.904.06	drill bit for quick coupling, Ø 2.5 x 125 mm, spiral length 60 mm, double spiral	1
2.906.40	pilot drill bit f.quick coupling, Ø 3.0 x 100 mm, bore length 60 mm, cannulated, I Ø 1.8 mm	1