



Calcaneus Plates


for Screws \varnothing 3.5 mm; Titanium and Stainless Steel



Item no. 5.862.01



Item no. 5.863.01



Item no. 4.862.01



Item no. 4.863.01

Item no.	Item description
----------	------------------

Calcaneus plates standard titanium

5.862.01	calcaneus plate 3.5, 10 holes, length 73, big
5.863.01	calcaneus plate 3.5, 8 holes, length 55, small

Calcaneus plates standard stainless steel

4.862.01	calcaneus plate 3.5, 10 holes, length 73, big
4.863.01	calcaneus plate 3.5, 8 holes, length 55, small