



# Variable Angle-stable Calcaneus Plates

for Screws  $\varnothing$  3.5 mm and  $\varnothing$  4.0 mm; Titanium



Item no. 5.8793.09



Item no. 5.895.00



Item no. 5.8813.01

Item no.

Item description

## Variable angle-stable calcaneus plates

5.8783.09

calcaneus plate with conic thread 3.5 and 4.0, 9 holes, length 65, right, titanium

5.8793.09

calcaneus plate with conic thread 3.5 and 4.0, 9 holes, length 65, left, titanium

5.894.00

calcaneus plate with conic thread 3.5 and 4.0, 10 holes, length 72, right, titanium

5.895.00

calcaneus plate with conic thread 3.5 and 4.0, 10 holes, length 72, left, titanium

5.8803.01

calcaneus plate with conic thread in 3 holes 3.5 and 4.0, 8 holes, angle-stable variable, right, titanium

5.8813.01

calcaneus plate with conic thread in 3 holes 3.5 and 4.0, 8 holes, angle-stable variable, left, titanium