

---

# Lateral Clavicle Plate 3.5 vas

Variable angle-stable for screws  $\varnothing$  3.5 mm; titanium

---



Item no. 5.3513.34, left



Item no. 5.3513.36, left

Item no.	Item description
----------	------------------

#### Variable angle-stable clavicle plate 3.5; titanium

5.3503.34	Clavicle plate with conical thread 3.5; 3/4 holes; 65 x 10 mm; lateral; vas; right; titanium
5.3503.36	Clavicle plate with conical thread 3.5; 3/6 holes; 82 x 10 mm; lateral; vas; right; titanium
5.3513.34	Clavicle plate with conical thread 3.5; 3/4 holes; 65 x 10 mm; lateral; vas; left; titanium
5.3513.36	Clavicle plate with conical thread 3.5; 3/6 holes; 82 x 10 mm; lateral; vas; left; titanium