



Fracture Compression Screws and Drill-Fix Forefoot Screws

for Forefoot Treatment / for Dens Screwing; Titanium



Fracture compression screw
Item no. 3.123.20



Drill-fix forefoot screw
Item no. 3.127.11

Item no.	Item description	EAN
Fracture compression screws Ø 3.0 mm / Ø 4.0 mm for forefoot treatment		
3.123.10	length 10 mm; thread 4; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088416
3.123.12	length 12 mm; thread 4; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088607
3.123.14	length 14 mm; thread 4; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088614
3.123.16	length 16 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088621
3.123.18	length 18 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088638
3.123.20	length 20 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088515
3.123.22	length 22 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088652
3.123.24	length 24 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088539
3.123.26	length 26 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088546
3.123.28	length 28 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088553
3.123.30	length 30 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088454
Fracture compression screws Ø 3.0 mm / Ø 4.0 mm for dens screwing		
3.123.32	length 32 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088805
3.123.34	length 34 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088669
3.123.36	length 36 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088676
3.123.38	length 38 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088683
3.123.40	length 40 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088690
3.123.42	length 42 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088706
3.123.44	length 44 mm; thread 6; self drilling; self tapping; cannulated; inside Ø 1.3 mm; titanium	4 047075 088560
Drill-fix forefoot screws Ø 2.0 mm		
3.127.11	length 11 mm; total length 26 mm; titanium	4 047075 088591
3.127.12	length 12 mm; total length 27 mm; titanium	4 047075 088836
3.127.13	length 13 mm; total length 28 mm; titanium	4 047075 088843
3.127.14	length 14 mm; total length 29 mm; titanium	4 047075 088850