

HTO plate 6.0/4.5 as medial

angle-stable; for screws Ø 6.0/4.5 mm; titanium



Item no. 5.7242.34



Item no. 5.7252.37



Item no. 5.7221.100



Item no. 3.721.09

Item no.	Item description	EAN
HTO plate 6.0/4.5 as medial		
5.7242.34	Königsee HTO plate medial with conic thread 6.0/4.5; 3/4 holes, 63 mm x 17 mm; right; titanium	4 047075 115907
5.7252.34	Königsee HTO plate medial with conic thread 6.0/4.5; 3/4 holes, 63 mm x 17 mm; left; titanium	4 047075 115914
5.7242.37	Königsee HTO plate medial with conic thread 6.0/4.5; 3/7 holes, 93 mm x 17 mm; right; titanium	4 047075 118045
5.7252.37	Königsee HTO plate medial with conic thread 6.0/4.5; 3/7 holes, 93 mm x 17 mm; left; titanium	4 047075 118038
HTO spacer		
5.7221.50	spacer for Königsee HTO plate, 8.0 mm x 5.0 mm, titanium	4 047075 117826
5.7221.75	spacer for Königsee HTO plate, 8.0 mm x 7.5 mm, titanium	4 047075 117918
5.7221.100	spacer for Königsee HTO plate, 8.0 mm x 10.0 mm, titanium	4 047075 117925
5.7221.125	spacer for Königsee HTO plate, 8.0 mm x 12.5 mm, titanium	4 047075 117932
5.7221.150	spacer for Königsee HTO plate, 8.0 mm x 15.0 mm, titanium	4 047075 117819
HTO screw for spacer		
3.721.09	fixing screw for distance holder, M4 x 9.0 mm, titanium	4 047075 110841