

## H-Plates Mini



H-plate Mini 1.5  
9 x 8.2 mm; 4 holes  
plate thickness 0.9 mm



H-plate Mini 2.0  
12 x 11 mm; 4 holes  
plate thickness 0.9 mm

## Straight plates Mini 2.0 as

Straight plate Mini 2.0 as  
plate thickness 1.2 mm; width 5 mm



5.9012.05: 5 holes; length 29 mm  
5.9012.13: 13 holes; length 77 mm

Straight plate Mini 2.0 as round bridge  
plate thickness 1.2 mm; 4 holes; 29 x 5 mm



## Contact details

If you are interested in the  
**HAND PLATES MINI**  
please feel free to contact us  
or your local sales representatives.



## Königsee Implantate GmbH

OT Aschau | Am Sand 4  
07426 Allendorf  
Germany

export@koenigsee-implantate.de  
+ 49 36738 498-560


All implants are available on request  
as non-sterile packed or sterile packed.

[www.koenigsee-implantate.de](http://www.koenigsee-implantate.de)

Copyright © Königsee Implantate GmbH/DaTi | Edition 01 | 2020-02

## HAND PLATES MINI

Titanium plates range



 **Königsee**  
Implantate

## T-Plates Mini 2.0 as

### T-Plate Mini 2.0 as narrow

plate thickness 1.2 mm; width 5 mm



5.9312.24: 2/4 holes; length 30 mm  
5.9312.26: 2/6 holes; length 42 mm  
5.9312.28: 2/8 holes; length 54 mm

### T-Plate Mini 2.0 as broad

grade; plate thickness 1.2 mm; width 5 mm



5.9332.34: 3/4 holes; length 30 mm  
5.9332.36: 3/6 holes; length 42 mm  
5.9332.38: 3/8 holes; length 54 mm

### T-Plate Mini 2.0 as Δ

3/4 holes; plate thickness 1.2 mm; 33 x 5 mm



### Y-Plate Mini 2.0 as

Y-form; plate thickness 1.2 mm; width 5 mm



5.9322.34: 3/4 holes; length 33 mm  
5.9322.36: 3/6 holes; length 45 mm  
5.9322.38: 3/8 holes; length 57 mm

## Cortical screw as

Ø 2.0 mm with depth stop; fully threaded; self tapping



 Internal hexagon

## Double bridge plate 2.0 as



5.9052.04S: 4 holes; 15 x 11 mm; oblique; plate thickness 1.2 mm

5.9052.07S: 7 holes; 25 x 10 mm; oblique; plate thickness 1.2 mm

5.9052.10: 10 holes; 45 x 11 mm; straight; plate thickness 1.2 mm

## K-Plate Mini 2.0 as

6 holes; plate thickness 1.2 mm; length 37 mm



## Straight plates Mini

### Straight plate Mini 1.5

20 holes; plate thickness 1.0 mm; 100 x 4 mm



### Straight plate Mini 2.0

with round holes; plate thickness 1.0 mm; width 5 mm



5.900.03: 3 holes; length 17 mm  
5.900.04: 4 holes; length 23 mm  
5.900.05: 5 holes; length 29 mm  
5.900.06: 6 holes; length 35 mm

## T-Plates Mini 2.0

### T-Plate Mini 2.0

plate thickness 1.5 mm; 50 x 5 mm



5.416.39: 3/9 holes  
5.416.49: 4/9 holes

### T-Plate Mini 2.0 long hole

long hole in shank; plate thickness 1.0 mm; width 5.5 mm



5.930.05: 5 holes; length 43 mm  
5.930.02: 2 holes; length 22 mm  
5.930.06: 6 holes; length 50 mm  
5.930.03: 3 holes; length 29 mm  
5.930.07: 7 holes; length 57 mm  
5.930.04: 4 holes; length 36 mm  
5.930.08: 8 holes; length 64 mm

### T-Plate Mini 2.0/2.7

with round holes 2.0 and 2.7; plate thickness 1.0 mm; width 7 mm



5.850.05: 5 holes; length 48 mm  
5.850.02: 2 holes; length 24 mm  
5.850.06: 6 holes; length 56 mm  
5.850.03: 3 holes; length 32 mm  
5.850.07: 7 holes; length 64 mm  
5.850.04: 4 holes; length 40 mm  
5.850.08: 8 holes; length 72 mm

## Cortical screw

Ø 2.0 mm; fully threaded; self tapping



 Internal hexagon