

Distal Humerus plate 3.5 as









radial/lateral with Extension; angle-stable; for screws \varnothing 3.5 mm; titanium



Item no. 5.8522.56



Item no. 5.8532.56

Item no.	Item description	EAN
5.8522.54	distal humerus plate with conical thread 3.5; 5/4 holes; length 78 mm; lateral with extension; right; titanium	 4 047075 160976
5.8522.56	distal humerus plate with conical thread 3.5; 5/6 holes; length 102 mm; lateral with extension; right; titanium	 4 047075 160983
5.8522.58	distal humerus plate with conical thread 3.5; 5/8 holes; length 126 mm; lateral with extension; right; titanium	 4 047075 160990
5.8522.510	distal humerus plate with conical thread 3.5; 5/10 holes; length 150 mm; lateral with extension; right; titanium	 4 047075 160969
5.8532.54	distal humerus plate with conical thread 3.5; 5/4 holes; length 78 mm; lateral with extension; left; titanium	 4 047075 161010
5.8532.56	distal humerus plate with conical thread 3.5; 5/6 holes; length 102 mm; lateral with extension; left; titanium	 4 047075 161027
5.8532.58	distal humerus plate with conical thread 3.5; 5/8 holes; length 126 mm; lateral with extension; left; titanium	 4 047075 161034
5.8532.510	distal humerus plate with conical thread 3.5; 5/10 holes; length 150 mm; lateral with extension; left; titanium	 4 047075 161003