

# Set Distal Humerus plate 3.5/4.5 as IH T

angle-stable; for screws Ø 3.5 mm/4.5 mm; titanium; Set no. 19.675.010




Item no. 5.8591.85



Item no. 5.8591.86





Item no. 5.8591.88







Set no. 19.675.010



Item no.	Item description	Quant.
<b>Implants – plates</b>		
<b>Distal humerus plates with conical thread 3.5 and 4.5; titanium</b>		
5.8591.85	Distal humerus plate with conic thread 3.5/4.5; 8/5 holes; length 141 mm	1
5.8591.86	Distal humerus plate with conic thread 3.5/4.5; 8/6 holes; length 169 mm	1
5.8591.88	Distal humerus plate with conic thread 3.5/4.5; 8/8 holes; length 199 mm	1

<b>Implants - screws</b>		
<b>Cortical screw with conical head thread; Ø 3.5 mm; fully threaded; self tapping; titanium</b>		
		
 Internal hexagon		
3.133.20	length 20 mm	5
3.133.22	length 22 mm	5
3.133.24	length 24 mm	5
3.133.26	length 26 mm	5
3.133.28	length 28 mm	5
3.133.30	length 30 mm	5
3.133.32	length 32 mm	5
3.133.34	length 34 mm	5
3.133.36	length 36 mm	5
3.133.38	length 38 mm	5
3.133.40	length 40 mm	5
3.133.42	length 42 mm	5
3.133.45	length 45 mm	5
3.133.48	length 48 mm	5
3.133.50	length 50 mm	5

Item no.	Item description	Quant.
<b>Implants - screws</b>		
<b>Cortical screw Ø 3.5 mm; fully threaded; self tapping; titanium</b>		
		
 Internal hexagon		
3.132.24	length 24 mm	2
3.132.26	length 26 mm	2
3.132.28	length 28 mm	2
3.132.30	length 30 mm	2
3.132.32	length 32 mm	2
3.132.34	length 34 mm	2
3.132.36	length 36 mm	2
3.132.38	length 38 mm	2
3.132.40	length 40 mm	2
<b>Cortical screw Ø 4.5 mm; fully threaded; self tapping; titanium</b>		
		
 Internal hexagon		
3.152.18	length 18 mm	2
3.152.20	length 20 mm	2
3.152.22	length 22 mm	2
3.152.24	length 24 mm	2
3.152.26	length 26 mm	2
3.152.28	length 28 mm	2
3.152.30	length 30 mm	2
3.152.32	length 32 mm	2
3.152.34	length 34 mm	2
3.152.36	length 36 mm	2
3.152.38	length 38 mm	2
3.152.40	length 40 mm	2
3.152.42	length 42 mm	2



Königsee Implantate GmbH  
OT Aschau | Am Sand 4  
07426 Allendorf | Germany

Fon +49 36738 498-560  
Fax +49 36738 498-569


Mail [export@koenigsee-implantate.de](mailto:export@koenigsee-implantate.de)

For further information please visit:  
[www.koenigsee-implantate.com](http://www.koenigsee-implantate.com)

HeBo | Edition 03 | 2022-10  
All pictured products are CE compliant.

# Set Distal Humerus plate 3.5/4.5 as IH T

angle-stable; for screws Ø 3.5 mm/4.5 mm; titanium; Set no. 19.675.010

Item no.	Item description	Quant.
<b>Implants - screws</b>		
	Cortical screw with conical head thread; Ø 4.5 mm; fully threaded; self tapping; titanium  Internal hexagon	
3.545.18	length 18 mm	3
3.545.20	length 20 mm	3
3.545.22	length 22 mm	3
3.545.24	length 24 mm	3
3.545.26	length 26 mm	3
3.545.28	length 28 mm	3
3.545.30	length 30 mm	3
3.545.32	length 32 mm	3
3.545.34	length 34 mm	3
3.545.36	length 36 mm	3
3.545.38	length 38 mm	3
3.545.40	length 40 mm	3
3.545.42	length 42 mm	3
3.545.45	length 45 mm	3
3.545.50	length 50 mm	3

Item no.	Item description	Quant.
<b>Instruments for small fragment</b>		
2.9406.25	hex screwdriver with handle for screws Ø 2.7mm - Ø 4.0 mm, AF 2.5 mm, conical	1
2.953.60	gauge with clasp for screws Ø 3.5mm, Ø 4.0mm, measuring range 60mm	1
2.904.06	drill bit for quick coupling, Ø 2.5 x 125 mm, spiral length 60 mm, double spiral	1
2.977.01	drill guide for angle stable screwing, length 45 mm, small fragment	2

<b>Instruments for large fragment</b>		
2.9406.35	hex screwdriver with handle for screws Ø 4.5 mm - Ø 7.0 mm, AF 3.5 mm, length 215 mm	1
2.953.90	gauge with clasp for screws Ø 4.5 mm, Ø 6.0 mm, measuring range 90 mm	1
2.904.09	drill bit for quick coupling, Ø 3.2 x 195 mm, spiral length 60 mm, double spiral	1
2.977.02	drill guide for angle stable screwing, length 70 mm, large fragment	2

<b>Instrument for large fragment and small fragment</b>		
2.954.01	screw forceps, self holding	1


<b>Wires</b>		
6.031.20	Kirschner wire with trocar point and round end, Ø 2.0 x 150 mm, stainless steel	5

<b>Container</b>		
19.674.010	perforated autoclavable container with inset for instruments/implants distal humerus plates, angle-stable, SF/LF, titanium	1

All implants are available on request as non-sterile packed or sterile packed.

# Set Distal Humerus plate 3.5/4.5 as IH T

angle-stable; for screws Ø 3.5 mm/4.5 mm; titanium; Set no. 19.675.010



Set no. 19.675.010



All implants are available on request as non-sterile packed or sterile packed.