

Set ImplantEx®

Set no. 19.511.00



Item no. 2.9366.15



Item no. 2.948.35A




Item no. 2.9101.55



Item no. 2.9042.050



Item no. 2.9481.45



◆ Internal hexagon
 ★ Sixlobe
 ✕ Cross head
 Set no. 19.511.00

Item no.	Item description	Quant.
Instruments – Screw driver		
Screwdriver with handle; Cross head		
2.9366.15	cross head; size 1.5	1

Screwdriver with handle; Sixlobe		
2.9386.06	SL 6	1
2.9386.07	SL 7	1
2.9386.08	SL 8	1
2.9386.10	SL 10	1
2.9386.15	SL 15	1
2.9386.25	SL 25	1

Screwdriver with handle; Internal hexagon		
2.9406.15	IH for screws Ø 1.5 – Ø 2.0; wrench size 1.5; kon.	1
2.9406.25	IH for screws Ø 2.7 – Ø 4.0; wrench size 2.5; kon.	1
2.9406.35	IH for screws Ø 4.5 – Ø 7.0; wrench size 3.5; L215	1

Instruments – Drill bits		
Drill bits; Round; DLC		
2.9042.025	2.5; Spiral 14; 2-flute; length 46	1
2.9042.030	3.0; Spiral 17; 2-flute; length 50	1
2.9042.040	4.0; Spiral 22; 2-flute; length 54	1
2.9042.050	5.0; Spiral 26; 2-flute; length 62	1

Container		
19.510.00	perforated autoclavable container with inset for IMPLANTEX®	1

Item no.	Item description	Quant.
Instruments – Removal tools		
Removal tools QC for IH		
2.948.15A	wrench size 1.5	1
2.948.25A	wrench size 2.5	1
2.948.35A	wrench size 3.5	1

Instruments – Burr hollows		
Burr hollows QC		
2.9101.15	for screws Ø 1.5 – Ø 2.0	1
2.9101.35	for screws Ø 2.7 – Ø 4.0	1
2.9101.45	for screws Ø 4.5 – Ø 7.0	1
2.9101.55	for screws Ø 7.0 – Ø 7.5	1

Instruments – Removal tools for screw shaft		
Removal tools QC for screw shaft		
2.9481.15	for screws Ø 1.5 – Ø 2.0	1
2.9481.24	for screws Ø 2.0 – Ø 2.7	1
2.9481.27	for screws Ø 2.7 – Ø 4.0	1
2.9481.45	for screws Ø 4.5 – Ø 6.0	1

Further instruments		
2.975.03	metal hammer; weight 300 g	1
10.2406.05	gouge	1
10.426.10	screw holding forceps; mini fragment	1
10.426.20	screw holding forceps; SF and LF	1
10.305.05	T-handle w. quick coupling; L 93.5	1



Königsee Implantate GmbH
 OT Aschau | Am Sand 4
 07426 Allendorf | Germany

Fon +49 36738 498-560
 Fax +49 36738 498-569

Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com


HeBo | Edition 02 | 2020-09
 All pictured products are CE compliant.

Set ImplantEx®


Set no. 19.511.00



Set no. 19.511.00 upper level



Set no. 19.511.00 lower level



◆ Internal hexagon
 ★ Sixlobe
 ✕ Cross head Set no. 19.511.00

Optional:

Screwdriver

- 2.9366.25 cross-head screwdriver with handle size 2.5
- 2.9366.35 cross-head screwdriver with handle size 3.5
- 10.266.02 hex screwdriver with handle;wrench width 2;length 350; cannulated;inside Ø1.3
- 2.9436.20 hex screwdriver w.handle f.fracture compression; cann.;|Ø1.3; wrench size 2;L215mm
- 2.9435.75 hex screwdriver with handle for Ø7.5;cannulated;|Ø2.9; wrench size 4;L 324mm

Removal tools QC for IH

- 2.948.50A screwdriver for extraction for quick coupling;wrench width 5.0
- 2.948.15 screwdriver with T-handle for extraction of screws with wrench size 1.5
- 2.948.25 screwdriver with T-handle for extraction of screws with wrench size 2.5
- 2.948.35 screwdriver with T-handle for extraction of screws with wrench size 3.5

Drill bits

- 2.9042.060 Drill bit round6.0;Spiral57;2-flute;L93;DLC
- 2.9042.070 Drill bit round7.0;Spiral69;2-flute;L109;DLC

Removal tools QC for screw shaft

- 2.9481.65 screwdriver for extraction from shank of screws Ø6.5 to Ø7.0 for quick coupling
- 2.9481.73 screwdriver for extraction from shank of screws Ø7.3 to Ø7.5 for quick coupling

Countersinks

- 10.505.15 countersink for dental coupling;length 57;for screws Ø1.5 and Ø2.0
- 10.505.40 countersink for quick coupling;length 72;for screws Ø2.7 / Ø3.5 / Ø4.0
- 10.505.81 countersink for quick coupling;length 110;for screws Ø4.5 and Ø6.5

Further instruments

- 10.421.06 extraction forceps for wires Ø1.0-2.0
- 10.421.07 extraction forceps for wires Ø2.0-3.0
- 10.305.10 coupling device from dental coupling to quick coupling



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498-560
Fax +49 36738 498-569

Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com

HeBo | Edition 02 | 2020-09
All pictured products are CE compliant.