

# KFN-Königsee Femoral Nail

Titanium



Item no.	Item description	EAN
----------	------------------	-----

## Implants - Nails

### Königsee Femoral Nail 125°; titanium

7.3KF.10/210	Königsee Femoral Nail 125°; Ø 10 x 210 mm; cannulated; inside Ø 4.7 mm	4 047075 131723
7.3KF.11/210	Königsee Femoral Nail 125°; Ø 11 x 210 mm; cannulated; inside Ø 4.7 mm	4 047075 131730
7.3KF.12/210	Königsee Femoral Nail 125°; Ø 12 x 210 mm; cannulated; inside Ø 4.7 mm	4 047075 131747
7.3KF.13/210	Königsee Femoral Nail 125°; Ø 13 x 210 mm; cannulated; inside Ø 4.7 mm	4 047075 131754
7.3KF10,5L/340	Königsee Femoral Nail 125°; Ø 10.5 x 340 mm; cannulated; inside Ø 4.7 mm; left	4 047075 142439
7.3KF10,5R/340	Königsee Femoral Nail 125°; Ø 10.5 x 340 mm; cannulated; inside Ø 4.7 mm; right	4 047075 142453
7.3KF10,5L/380	Königsee Femoral Nail 125°; Ø 10.5 x 380 mm; cannulated; inside Ø 4.7 mm; left	4 047075 142446
7.3KF10,5R/380	Königsee Femoral Nail 125°; Ø 10.5 x 380 mm; cannulated; inside Ø 4.7 mm; right	4 047075 142460

### Königsee Femoral Nail 130°; titanium

7.1KF.10/210	Königsee Femoral Nail 130°; Ø 10 x 210 mm; cannulated; inside Ø 4.7 mm	4 047075 129140
7.1KF.11/210	Königsee Femoral Nail 130°; Ø 11 x 210 mm; cannulated; inside Ø 4.7 mm	4 047075 129157
7.1KF.12/210	Königsee Femoral Nail 130°; Ø 12 x 210 mm; cannulated; inside Ø 4.7 mm	4 047075 129164
7.1KF.13/210	Königsee Femoral Nail 130°; Ø 13 x 210 mm; cannulated; inside Ø 4.7 mm	4 047075 135165
7.1KF10,5L/340	Königsee Femoral Nail 130°; Ø 10.5 x 340 mm; cannulated; inside Ø 4.7 mm; left	4 047075 129171
7.1KF10,5R/340	Königsee Femoral Nail 130°; Ø 10.5 x 340 mm; cannulated; inside Ø 4.7 mm; right	4 047075 129195
7.1KF10,5L/420	Königsee Femoral Nail 130°; Ø 10.5 x 420 mm; cannulated; inside Ø 4.7 mm; left	4 047075 129188
7.1KF10,5R/420	Königsee Femoral Nail 130°; Ø 10.5 x 420 mm; cannulated; inside Ø 4.7 mm; right	4 047075 129201

## Implants - Screws; titanium alloy; TiOB®-coated; bioinert

3.157.25	interlocking screw for KFN; Ø 5.5 x 25 mm; self tapping	4 047075 093199
3.157.26	interlocking screw for KFN; Ø 5.5 x 26 mm; self tapping	4 047075 049530
3.157.28	interlocking screw for KFN; Ø 5.5 x 28 mm; self tapping	4 047075 062454
3.157.30	interlocking screw for KFN; Ø 5.5 x 30 mm; self tapping	4 047075 093083
3.157.32	interlocking screw for KFN; Ø 5.5 x 32 mm; self tapping	4 047075 093090



Königsee Implantate GmbH  
 OT Aschau | Am Sand 4  
 07426 Allendorf | Germany

Fon +49 36738 498-560  
 Fax +49 36738 498-569

Mail [export@koenigsee-implantate.de](mailto:export@koenigsee-implantate.de)

For further information please visit:  
[www.koenigsee-implantate.com](http://www.koenigsee-implantate.com)

HeBo | Edition 02 | 2020-02  
 All pictured products are CE compliant.

# KFN-Königsee Femoral Nail

Titanium



Item no. 7.1KF.10/210



Item no. 7.1KF.11/210



Item no. 7.1KF.12/210



3.157.34



3.428.110



7.311.12



7.311.30



7.311.01



3.310.52



7.313.120

Item no.	Item description	EAN
----------	------------------	-----

## Implants - Screws; titanium alloy; TiOB<sup>®</sup>-coated; bioinert

3.157.34	interlocking screw for KFN; Ø 5.5 x 34 mm; self tapping	4 047075 093106
3.157.35	interlocking screw for KFN; Ø 5.5 x 35 mm; self tapping	4 047075 093113
3.157.36	interlocking screw for KFN; Ø 5.5 x 36 mm; self tapping	4 047075 093120
3.157.38	interlocking screw for KFN; Ø 5.5 x 38 mm; self tapping	4 047075 093014
3.157.40	interlocking screw for KFN; Ø 5.5 x 40 mm; self tapping	4 047075 093021
3.157.42	interlocking screw for KFN; Ø 5.5 x 42 mm; self tapping	4 047075 093144
3.157.44	interlocking screw for KFN; Ø 5.5 x 44 mm; self tapping	4 047075 093045
3.157.45	interlocking screw for KFN; Ø 5.5 x 45 mm; self tapping	4 047075 093052
3.157.50	interlocking screw for KFN; Ø 5.5 x 50 mm; self tapping	4 047075 148523
3.157.55	interlocking screw for KFN; Ø 5.5 x 55 mm; self tapping	4 047075 148530
3.157.60	interlocking screw for KFN; Ø 5.5 x 60 mm; self tapping	4 047075 148547
3.157.65	interlocking screw for KFN; Ø 5.5 x 65 mm; self tapping	4 047075 148554
3.428.65	anti-rotation screw for KFN; Ø 6.0 x 65 mm	4 047075 121021
3.428.70	anti-rotation screw for KFN; Ø 6.0 x 70 mm	4 047075 121038
3.428.75	anti-rotation screw for KFN; Ø 6.0 x 75 mm	4 047075 097029
3.428.80	anti-rotation screw for KFN; Ø 6.0 x 80 mm	4 047075 097036
3.428.85	anti-rotation screw for KFN; Ø 6.0 x 85 mm	4 047075 097180
3.428.90	anti-rotation screw for KFN; Ø 6.0 x 90 mm	4 047075 097197
3.428.95	anti-rotation screw for KFN; Ø 6.0 x 95 mm	4 047075 097050
3.428.100	anti-rotation screw for KFN; Ø 6.0 x 100 mm	4 047075 097135
3.428.105	anti-rotation screw for KFN; Ø 6.0 x 105 mm	4 047075 097142
3.428.110	anti-rotation screw for KFN; Ø 6.0 x 110 mm	4 047075 097159
3.428.115	anti-rotation screw for KFN; Ø 6.0 x 115 mm	4 047075 097005
3.428.120	anti-rotation screw for KFN; Ø 6.0 x 120 mm	4 047075 097012

# KFN-Königsee Femoral Nail

Titanium



Item no. 7.1KF.10/210



Item no. 7.1KF.11/210



Item no. 7.1KF.12/210



3.157.34



3.428.110



7.311.12



7.311.30



7.311.01



3.310.52



7.313.120

Item no.	Item description	EAN
----------	------------------	-----

## Implants - Screws

7.311.12	locking screw without arresting bolt for KFN; length 28 mm; M 12; titanium alloy; TiOB <sup>®</sup> -coated; bioinert	
7.311.30	locking screw with arresting bolt for KFN; length 30 mm; M 12; titanium alloy; TiOB <sup>®</sup> -coated; bioinert	
7.311.01	compression screw for DHS; M 4 x 27 mm; length 35 mm; SW 3.5; titanium alloy	
3.310.52	supporting screw sleeve for KFN; Ø 12 x 52 mm; titanium alloy; TiOB <sup>®</sup> -coated; bioinert	

## Implants - Screws; titanium alloy; TiOB<sup>®</sup>-coated; bioinert

7.313.70	supporting screw for KFN; Ø 10 x 70 mm	
7.313.75	supporting screw for KFN; Ø 10 x 75 mm	
7.313.80	supporting screw for KFN; Ø 10 x 80 mm	
7.313.85	supporting screw for KFN; Ø 10 x 85 mm	
7.313.90	supporting screw for KFN; Ø 10 x 90 mm	
7.313.95	supporting screw for KFN; Ø 10 x 95 mm	
7.313.100	supporting screw for KFN; Ø 10 x 100 mm	
7.313.105	supporting screw for KFN; Ø 10 x 105 mm	
7.313.110	supporting screw for KFN; Ø 10 x 110 mm	
7.313.115	supporting screw for KFN; Ø 10 x 115 mm	
7.313.120	supporting screw for KFN; Ø 10 x 120 mm	
7.313.125	supporting screw for KFN; Ø 10 x 125 mm	
7.313.130	supporting screw for KFN; Ø 10 x 130 mm	

TiOB<sup>®</sup> is a trademark of Königsee Implantate GmbH registered in the EU countries, Switzerland and the USA.

It is a surface coating that brings optimized properties.



Königsee Implantate GmbH  
OT Aschau | Am Sand 4  
07426 Allendorf | Germany

Fon +49 36738 498-560  
Fax +49 36738 498-569

Mail [export@koenigsee-implantate.de](mailto:export@koenigsee-implantate.de)

For further information please visit:  
[www.koenigsee-implantate.com](http://www.koenigsee-implantate.com)

HeBo | Edition 02 | 2020-02  
All pictured products are CE compliant.