



Set KFN-Königsee Femoral Nail

Titanium; Set no. 19.344.00





Item no. 7.1KF.10/210 Item no. 7.1KF.11/210 Item no. 7.1KF.12/210





Container Implants; Set no. 19.345.00  Internal hexagon


Item no.	Item description	Quant.
Implants – KFN nails		
KFN-Königsee Femoral Nail 130°; titanium		
7.1KF.10/210	Ø10x210mm; cannulated; inside-Ø4.7 mm	1
7.1KF.11/210	Ø11x210mm; cannulated; inside-Ø4.7 mm	1
7.1KF.12/210	Ø12x210mm; cannulated; inside-Ø4.7 mm	1


Implants – screws		
Interlocking screw for KFN Ø 5.5 mm; self tapping; titanium alloy; TiOB®-coated; bioinert 		
3.157.26	length 26 mm	2
3.157.28	length 28 mm	2
3.157.30	length 30 mm	2
3.157.32	length 32 mm	2
3.157.34	length 34 mm	2
3.157.36	length 36 mm	2
3.157.38	length 38 mm	2
3.157.40	length 40 mm	2
3.157.42	length 42 mm	2
3.157.44	length 44 mm	2

Locking screw for KFN; titanium alloy; TiOB®-coated; bioinert 		
7.311.12	without arresting bolt; length 28 mm; M 12	2
7.311.30	with arresting bolt; length 30 mm; M 12	2

Compression screw for KFN; titanium alloy 		
7.311.01	thread 27; M 4; L 35 mm; wrench width 3.5	2

Sleeve of support. screw; tit. alloy; TiOB®-coated; bioinert 		
3.310.52	supporting screw sleeve f. KFN; Ø 12 x 52 mm	2


Item no.	Item description	Quant.
Implants – screws		
Anti-rotation screw for KFN; Ø 6.0 mm; titanium alloy; TiOB®-coated; bioinert 		
3.428.75	length 75 mm	1
3.428.80	length 80 mm	1
3.428.85	length 85 mm	1
3.428.90	length 90 mm	1
3.428.95	length 95 mm	1
3.428.100	length 100 mm	1
3.428.105	length 105 mm	1
3.428.110	length 110 mm	1
3.428.115	length 115 mm	1
3.428.120	length 120 mm	1

Supporting screw for KFN; Ø 10 mm; titanium alloy; TiOB®-coated; bioinert 		
7.313.85	length 85 mm	1
7.313.90	length 90 mm	1
7.313.95	length 95 mm	1
7.313.100	length 100 mm	1
7.313.105	length 105 mm	1
7.313.110	length 110 mm	1
7.313.115	length 115 mm	1
7.313.120	length 120 mm	1
7.313.125	length 125 mm	1
7.313.130	length 130 mm	1

Container		
19.345.00	perforated autoclavable container with inset for implants KFN	1

Set KFN-Königsee Femoral Nail

Titanium; Set no. 19.344.00



Container 1 Instruments; Set no. 19.347.00  Internal hexagon

Item no.	Item description	Quant.
Instruments – color code red		
10.264.19	intramedullary reamer; Ø 17 x 269 mm; cannulated; inside Ø 3.5 mm	1
10.264.20	tissue protecting drill sleeve for intramedullary reamer; Ø 20 x 107 mm; inside Ø 17.2 mm	1
10.264.21	sleeve for intramedullary reamer; Ø 17 mm/ Ø 3.2 x 121 mm	1
10.264.22	trocars for sleeve for intramedullary reamer; Ø 3.0 x 132 mm	1

Wire – color code red		
2.984.30/400	guide wire; Ø 3.0 x 400 mm; thread 10; stainless steel	5

Item no.	Item description	Quant.
Instruments – color code yellow		
10.264.01	screwdriver for insertion	1
10.264.02	fixation rod for screwdriver for insertion; Ø 4.0 x 336 mm	1
10.264.03	hex screwdriver with handle and holding sleeve; wrench width 3.5	1
10.264.04	gauge for supporting screw	1
10.264.05	sleeve for supporting screw; Ø 16 mm/ Ø 14 x 200 mm	1
10.264.06	sleeve for supporting screw; Ø 14 mm/ Ø 2.7 x 211 mm	1
10.264.07	trocars for sleeve for supporting screw; Ø 2.5 x 222 mm	1
10.264.08	step drill-set for supporting screw; complete	1
10.264.27	T-handle with quick coupling for KFN; cannulated; inside Ø 4.2 mm	1

Wire – color code yellow		
2.984.25/400	guide wire; Ø 2.5 x 400 mm; thread 10; stainless steel	5

Container		
19.347.00	perforated autoclavable container with inset for instruments 1 KFN	1

TiOB® is a trademark of Königsee Implantate GmbH registered in the EU countries, Switzerland and the USA.

It is a surface treatment that brings optimized properties.



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498-560
Fax +49 36738 498-569


Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com

HeBo | Edition 02 | 2020-02
All pictured products are CE compliant.

Set KFN-Königsee Femoral Nail

Titanium; Set no. 19.344.00



Container 2 Instruments; Set no. 19.348.00  Internal hexagon

Item no.	Item description	Quant.
Instruments – color code green		
10.264.12	sleeve for anti-rotation screw; Ø 8.8 mm/ Ø 6.7 x 198 mm	1
10.264.14	sleeve for anti-rotation screw; Ø 6.5 mm/ Ø 2.2 x 205 mm	1
10.264.15	trocax for sleeve for anti-rotation screw; Ø 2.0 x 215 mm	1
10.264.16	drill bit; Ø 6.0 x 400 mm; quadruple spiral; spiral length 55 mm; cannulated; inside Ø 2.2 mm	1
10.264.17	hex screwdriver with handle; cannulated; inside Ø 2.2 mm; wrench width 3.5 mm	1
10.264.18	unscrewable sleeve; Ø 10 x 200.5 mm	1

Wire – color code green		
2.984.20/400	guide wire; Ø 2.0 x 400 mm; thread 10; stainless steel	5

Item no.	Item description	Quant.
Instruments – color code silver		
10.261.011	drill guide; Ø 8.0 x 150 mm; indented; silver	1
10.261.031	drill guide; Ø 4.8 x 160 mm; silver	1
10.261.15	drill bit; Ø 4.8 x 305 mm; spiral length 60 mm	1
10.261.16	gauge with clasp; measuring range 110 mm	1
2.930.35	hex screwdriver with T-handle for screw Ø 4.5 mm until Ø 7.0 mm; length 232 mm	1
10.263.01	impactor and extractor for KFN	1

Container		
19.348.00	perforated autoclavable container with inset for instruments 2 KFN	1

TiOB® is a trademark of Königsee Implantate GmbH registered in the EU countries, Switzerland and the USA.

It is a surface treatment that brings optimized properties.



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498-560
Fax +49 36738 498-569


Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com

HeBo | Edition 02 | 2020-02
All pictured products are CE compliant.

Set KFN-Königsee Femoral Nail

Titanium; Set no. 19.344.00



Container 3 Instruments; Set no. 19.349.00  Internal hexagon

Item no.	Item description	Quant.
Instruments – guide frame system		
10.264.50	guide frame KFN	1
10.264.50/3	locking screw for guide frame KFN	1
10.264.51	cap for guide frame; distal; complete	1
10.264.52	cap for guide frame 130°; complete	1
10.264.55	cap for guide frame 125°; complete	1

Item no.	Item description	Quant.
Container		
19.349.00	perforated autoclavable container with inset for instruments 3 KFN	1

Instruments		
10.262.04	bent awl; big	1
10.263.05	open-end wrench with handle; wrench width 17	1
10.264.09	screwdriver for extraction with t-handle	1
10.264.10	fixation rod for screwdriver for extraction; Ø 4.0 x 270 mm	1
10.264.11A	adapter for impactor for guide frame	1
10.264.23	tissue protector	1
10.264.24	plastic hammer	1
10.264.29	hex screwdriver with T-handle; wrench width 8	1
10.264.30/3	locking screw for guide frame KFN	1

TiOB® is a trademark of Königsee Implantate GmbH registered in the EU countries, Switzerland and the USA.

It is a surface treatment that brings optimized properties.



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498-560
Fax +49 36738 498-569


Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com


HeBo | Edition 02 | 2020-02
All pictured products are CE compliant.

Set KFN-Königsee Femoral Nail

Titanium; Set no. 19.344.00



Container implants KFN; set no. 19.345.00  Internal hexagon



Container 1 instruments; set no. 19.347.00  Internal hexagon

Set KFN-Königsee Femoral Nail

Titanium; Set no. 19.344.00



Container 2 instruments; set no. 19.348.00  Internal hexagon



Container 3 instruments; set no. 19.349.00  Internal hexagon