

Set LC titanium strip T

with clamp; 2 levels; titanium


not sterile packed Set no. 19.519.00



Item no. 10.271.20
Item no. 7.440.01
Item no. S7.440.01



Item no. 10.271.20/01




Set no. 19.519.00
not sterile packed

All implants mentioned are available individually, not sterile and sterile packed.

Item no.	Item no. sterile	Item description	Quant.
Implants			
LC titanium strip			
7.440.01	S7.440.01	LC titanium strip with clamp; 8.0 mm x 0.4 mm; length 260 mm	4
Instruments			
10.271.11/65	–	Strip loop instrument; 1/4"-coupling; straight; Ø 65 mm; length 225 mm	1
10.271.12/65	–	Strip loop instrument; 1/4"-coupling; beveled; Ø 65 mm; length 185 mm	1
10.271.11/45	–	Strip loop instrument; 1/4"-coupling; straight; Ø 45 mm; length 205 mm	1
10.271.12/45	–	Strip loop instrument; 1/4"-coupling; beveled; Ø 45 mm; length 175 mm	1
10.307.03	–	Handle with 1/4" coupling; length 162 mm	1
10.271.20	–	Universal clamping device for strips up to 8.0 mm x 0.4 mm	1
10.271.03	–	Pestle for LC titanium strip 8.0 mm x 0.4 mm; 1/4"-coupling; length 125 mm	1
10.271.05	–	Power cutter for LC titanium strip 8.0 mm x 0.4 mm; length 260 mm	1
2.975.03	–	Metal hammer; weight 300 g	1
Container			
19.519.10	–	Perforated autoclavable container with inset for instruments for titanium strips	1
Optional			
10.271.20/01	–	Cap for spindle handle for universal clamping device for strips up to 8.0 mm x 0.4 mm	1

Set LC titanium strip T
with clamp; 2 levels; titanium
not sterile packed Set no. 19.519.00



Set no. 19.519.00
not sterile packed