

# Set LC-Titanium Strip

with clamp; 2 levels; titanium; set-no. 19.519.00



Item no. 10.271.20



Set no. 19.519.00

Item no.	Item description	Quant.
<b>Implants</b>		
<b>LC-titanium strip with clamp</b>		
7.440.01	8.0mmx0.4mm; L 260mm	4

Item no.	Item description	Quant.
<b>Container</b>		
19.519.10	perforated autocl.container w.inset f.instruments for titanium strips	1


<b>Instruments</b>		
10.271.11/65	strip loop instrument; 1/4"-coupling; straight; Ø 65mm; L 225mm	1
10.271.12/65	strip loop instrument; 1/4"-coupling; beveled; Ø 65mm; L 185mm	1
10.271.11/45	strip loop instrument; 1/4"-coupling; straight; Ø 45mm; L 205mm	1
10.271.12/45	strip loop instrument; 1/4"-coupling; beveled; Ø 45mm; L 175mm	1
10.307.03	Handle with 1/4" coupling,length 162	1
10.271.20	universal clamping device f.strips up to 8mmx0.4mm/slopes up to Ø2mm; L220mm	1
10.271.03	pestle f.LC titanium strip 8.0mmx0.4mm; 1/4"-coupling; L 125mm	1
10.271.05	power cutter for LC titanium strip 8.0mmx0.4mm; L 260mm	1
2.975.03	metal hammer,weight 300 g	1

---

# Set LC-Titanium Strip

with clamp; 2 levels; titanium; set-no. 19.519.00

---



Set no. 19.519.00  
2 Level view