

Low contact plate 4.5/6.5

for screws Ø 4.5 mm and Ø 6.5 mm; stainless steel



Item no. 4.230.06



Item no. 4.240.08



Item no. 4.230.12; back

Item no.	Item description	EAN
----------	------------------	-----

Low contact plate 4.5/6.5 narrow



4.230.02	Low-contact plate 4.5 and 6.5; 2 holes; 34 x 12 mm; narrow; stainless steel	
4.230.03	Low-contact plate 4.5 and 6.5; 3 holes; 52 x 12 mm; narrow; stainless steel	
4.230.04	Low-contact plate 4.5 and 6.5; 4 holes; 70 x 12 mm; narrow; stainless steel	
4.230.05	Low-contact plate 4.5 and 6.5; 5 holes; 88 x 12 mm; narrow; stainless steel	
4.230.06	Low-contact plate 4.5 and 6.5; 6 holes; 106 x 12 mm; narrow; stainless steel	
4.230.07	Low-contact plate 4.5 and 6.5; 7 holes; 124 x 12 mm; narrow; stainless steel	
4.230.08	Low-contact plate 4.5 and 6.5; 8 holes; 142 x 12 mm; narrow; stainless steel	
4.230.09	Low-contact plate 4.5 and 6.5; 9 holes; 160 x 12 mm; narrow; stainless steel	
4.230.10	Low-contact plate 4.5 and 6.5; 10 holes; 178 x 12 mm; narrow; stainless steel	
4.230.11	Low-contact plate 4.5 and 6.5; 11 holes; 196 x 12 mm; narrow; stainless steel	
4.230.12	Low-contact plate 4.5 and 6.5; 12 holes; 214 x 12 mm; narrow; stainless steel	
4.230.13	Low-contact plate 4.5 and 6.5; 13 holes; 232 x 12 mm; narrow; stainless steel	
4.230.14	Low-contact plate 4.5 and 6.5; 14 holes; 250 x 12 mm; narrow; stainless steel	
4.230.15	Low-contact plate 4.5 and 6.5; 15 holes; 268 x 12 mm; narrow; stainless steel	
4.230.16	Low-contact plate 4.5 and 6.5; 16 holes; 286 x 12 mm; narrow; stainless steel	

Low contact plate 4.5/6.5 broad

4.240.06	Low-contact plate 4.5 and 6.5; 6 holes; 106 x 16 mm; broad; staggered holes; stainless steel	
4.240.07	Low-contact plate 4.5 and 6.5; 7 holes; 124 x 16 mm; broad; staggered holes; stainless steel	
4.240.08	Low-contact plate 4.5 and 6.5; 8 holes; 142 x 16 mm; broad; staggered holes; stainless steel	
4.240.09	Low-contact plate 4.5 and 6.5; 9 holes; 160 x 16 mm; broad; staggered holes; stainless steel	
4.240.10	Low-contact plate 4.5 and 6.5; 10 holes; 178 x 16 mm; broad; staggered holes; stainless steel	
4.240.11	Low-contact plate 4.5 and 6.5; 11 holes; 196 x 16 mm; broad; staggered holes; stainless steel	

Low contact plate 4.5/6.5

for screws Ø 4.5 mm and Ø 6.5 mm; stainless steel

Item no.	Item description	EAN
4.240.12	Low-contact plate 4.5 and 6.5; 12 holes; 214 x 16 mm; broad; staggered holes; stainless steel	 4 047075 103287
4.240.13	Low-contact plate 4.5 and 6.5; 13 holes; 232 x 16 mm; broad; staggered holes; stainless steel	 4 047075 103416
4.240.14	Low-contact plate 4.5 and 6.5; 14 holes; 250 x 16 mm; broad; staggered holes; stainless steel	 4 047075 103300
4.240.15	Low-contact plate 4.5 and 6.5; 15 holes; 268 x 16 mm; broad; staggered holes; stainless steel	 4 047075 103317
4.240.16	Low-contact plate 4.5 and 6.5; 16 holes; 286 x 16 mm; broad; staggered holes; stainless steel	 4 047075 103324
4.240.17	Low-contact plate 4.5 and 6.5; 17 holes; 304 x 16 mm; broad; staggered holes; stainless steel	 4 047075 103218
4.240.18	Low-contact plate 4.5 and 6.5; 18 holes; 322 x 16 mm; broad; staggered holes; stainless steel	 4 047075 103409
4.240.20	Low-contact plate 4.5 and 6.5; 20 holes; 358 x 16 mm; broad; staggered holes; stainless steel	 4 047075 103423