

# Mini plates hand 2.0 ws

angle-stable; for screws  $\varnothing$  2.0 mm; plate thickness 1.2 mm; titanium  
non-sterile packaging / sterile packaging



Item no. 5.9012.05  
Item no. S5.9012.05



Item no. 5.9022.04  
Item no. S5.9022.04



Item no. 5.9052.04S  
Item no. S5.9052.04S



Item no. 5.9052.07S  
Item no. S5.9052.07S



Item no. 5.9012.13  
Item no. S5.9012.13



Item no. 5.9032.06  
Item no. S5.9032.06



Item no. 5.9052.10  
Item no. S5.9052.10



Item no. 5.9312.26  
Item no. S5.9312.26

Item no.	EAN	Item description	Item no. sterile	EAN sterile
<b>Straight plate Mini 2.0 as</b>				
5.9012.05		Straight plate with conical thread $\varnothing$ 2.0 mm; 5 holes; 29 mm x 5.0 mm; titanium	S5.9012.05	
5.9012.13		Straight plate with conical thread $\varnothing$ 2.0 mm; 13 holes; 77 mm x 5.0 mm; titanium	S5.9012.13	
<b>Straight plate Mini 2.0 as bridge</b>				
5.9022.04		Mini bridge plate with conical thread $\varnothing$ 2.0 mm; 4 holes; 29 mm x 5.0 mm; titanium	S5.9022.04	
<b>K plate Mini 2.0 as</b>				
5.9032.06		K-plate with conical thread $\varnothing$ 2.0 mm; 6 holes; length 37 mm; titanium	S5.9032.06	
<b>Double bridge plate 2.0 as skew</b>				
5.9052.04S		Double bridge plate with conical thread $\varnothing$ 2.0 mm; 4 holes; 15 mm x 11 mm; oblique; titanium	S5.9052.04S	
5.9052.07S		Double bridge plate with conical thread $\varnothing$ 2.0 mm; 7 holes; 25 mm x 10 mm; oblique; titanium	S5.9052.07S	
5.9052.10		Double bridge plate with conical thread $\varnothing$ 2.0 mm; 10 holes; 45 mm x 11 mm; straight; titanium	S5.9052.10	
<b>T plate Mini 2.0 as narrow</b>				
5.9312.24		Mini T-plate with conical thread $\varnothing$ 2.0/2.0 mm; 2/4 holes; 30 mm x 5.0 mm; titanium	S5.9312.24	
5.9312.26		Mini T-plate with conical thread $\varnothing$ 2.0/2.0 mm; 2/6 holes; 42 mm x 5.0 mm; titanium	S5.9312.26	
5.9312.28		Mini T-plate with conical thread $\varnothing$ 2.0/2.0 mm; 2/8 holes; 54 mm x 5.0 mm; titanium	S5.9312.28	

# Mini plates hand 2.0 ws

angle-stable; for screws  $\varnothing$  2.0 mm; plate thickness 1.2 mm; titanium  
non-sterile packaging / sterile packaging



Item no. 5.9332.36  
Item no. S5.9332.36



Item no. 5.9322.36  
Item no. S5.9322.36



Item no. 5.9342.34  
Item no. S5.9342.34

Item no.	EAN	Item description	Item no. sterile	EAN sterile
<b>T plate Mini 2.0 as broad</b>				
5.9332.34		Mini T-plate with conical thread $\varnothing$ 2.0 mm/2.0 mm; 3/4 holes; 30 mm x 5.0 mm; straight; titanium	S5.9332.34	
5.9332.36		Mini T-plate with conical thread $\varnothing$ 2.0 mm/2.0 mm; 3/6 holes; 42 mm x 5.0 mm; straight; titanium	S5.9332.36	
5.9332.38		Mini T-plate with conical thread $\varnothing$ 2.0 mm/2.0 mm; 3/8 holes; 54 mm x 5.0 mm; straight; titanium	S5.9332.38	
<b>Y-plate Mini 2.0 as</b>				
5.9322.34		Mini T-plate with conical thread $\varnothing$ 2.0 mm/2.0 mm; 3/4 holes; 33 mm x 5.0 mm; y-form; titanium	S5.9322.34	
5.9322.36		Mini T-plate with conical thread $\varnothing$ 2.0 mm/2.0 mm; 3/6 holes; 45 mm x 5.0 mm; y-form; titanium	S5.9322.36	
5.9322.38		Mini T-plate with conical thread $\varnothing$ 2.0 mm/2.0 mm; 3/8 holes; 57 mm x 5.0 mm; y-form; titanium	S5.9322.38	
<b>T plate Mini 2.0 as</b>				
5.9342.34		Mini T-plate with conical thread 2.0 mm/2.0 mm; 3/4 holes; 33 mm x 5.0 mm; titanium	S5.9342.34	