

Set MF plate as IH T SET / IMPLANTS sterile

angle-stable; for screws Ø 2.0 mm; titanium

non-sterile packaging Set no. 19.508.00 / sterile packaging Set no. 19.508.00PS+



Item no. 5.9012.05
Item no. S5.9012.05



Item no. 5.9012.13
Item no. S5.9012.13



Item no. 5.9052.10
Item no. S5.9052.10



Item no. 5.9052.04S
Item no. S5.9052.04S



Item no. 5.9312.24
Item no. S5.9312.24



Item no. 5.9332.34
Item no. S5.9332.34



Item no. 5.9322.34
Item no. S5.9322.34



All listed implants are available individually in non-sterile and sterile packaging.

Implants in sets, supplied **non-sterile** Set no. 19.508.00

Implants in sets, supplied **sterile** Set no. 19.508.00PS+

Internal hexagon
Set no. 19.508.00
non-sterile packaging

Item no.	Item no. sterile	Item description	Quant.
Implantats – plates			
Straight plate Mini 2.0 as			
5.9012.05	S5.9012.05	Straight plate with conical thread Ø 2.0 mm; 5 holes; 29 mm x 5.0 mm; titanium	1
5.9012.13	S5.9012.13	Straight plate with conical thread Ø 2.0 mm; 13 holes; 77 mm x 5.0 mm; titanium	1
Double bridge plate 2.0 as			
5.9052.10	S5.9052.10	Double bridge plate with conical thread Ø 2.0 mm; 10 holes; 45 mm x 11 mm; straight; titanium	1
5.9052.04S	S5.9052.04S	Double bridge plate with conical thread Ø 2.0 mm; 4 holes; 15 mm x 11 mm; oblique; titanium	1
T plate Mini 2.0 as narrow			
5.9312.24	S5.9312.24	Mini T-plate with conical thread Ø 2.0/2.0; 2/4 holes; 30 mm x 5.0 mm; titanium	1
T plate Mini 2.0 as broad			
5.9332.34	S5.9332.34	Mini T-plate with conical thread Ø 2.0/2.0; 3/4 holes; 30 mm x 5.0 mm; straight; titanium	1
Y-plate Mini 2.0 as			
5.9322.34	S5.9322.34	Mini T-plate with conical thread Ø 2.0/2.0; 3/4 holes; 33 mm x 50 mm; y-form; titanium	1
Implants – screws			
Cortical screw 2.0 as st			
3.114.08	S3.114.08	Cortical screw; Ø 2.0 mm x 8.0 mm; angle-stable; fully threaded; self tapping; internal hexagon 1.5 mm; titanium	5
3.114.09	S3.114.09	Cortical screw; Ø 2.0 mm x 9.0 mm; angle-stable; fully threaded; self tapping; internal hexagon 1.5 mm; titanium	5

Set MF plate as IH T SET / IMPLANTS sterile

angle-stable; for screws Ø 2.0 mm; titanium

non-sterile packaging Set no. 19.508.00 / sterile packaging Set no. 19.508.00PS+

Item no.	Item no. sterile	Item description	Quant.
Implants – screws			
3.114.10	S3.114.10	Cortical screw; Ø 2.0 mm x 10 mm; angle-stable; fully threaded; self tapping; internal hexagon 1.5 mm; titanium	5
3.114.11	S3.114.11	Cortical screw; Ø 2.0 mm x 11 mm; angle-stable; fully threaded; self tapping; internal hexagon 1.5 mm; titanium	5
3.114.12	S3.114.12	Cortical screw; Ø 2.0 mm x 12 mm; angle-stable; fully threaded; self tapping; internal hexagon 1.5 mm; titanium	5
3.114.13	S3.114.13	Cortical screw; Ø 2.0 mm x 13 mm; angle-stable; fully threaded; self tapping; internal hexagon 1.5 mm; titanium	5
3.114.14	S3.114.14	Cortical screw; Ø 2.0 mm x 14 mm; angle-stable; fully threaded; self tapping; internal hexagon 1.5 mm; titanium	5
3.114.15	S3.114.15	Cortical screw; Ø 2.0 mm x 15 mm; angle-stable; fully threaded; self tapping; internal hexagon 1.5 mm; titanium	5
3.114.16	S3.114.16	Cortical screw; Ø 2.0 mm x 16 mm; angle-stable; fully threaded; self tapping; internal hexagon 1.5 mm; titanium	5
3.114.18	S3.114.18	Cortical screw; Ø 2.0 mm x 18 mm; angle-stable; fully threaded; self tapping; internal hexagon 1.5 mm; titanium	5
3.114.20	S3.114.20	Cortical screw; Ø 2.0 mm x 20 mm; angle-stable; fully threaded; self tapping; internal hexagon 1.5 mm; titanium	5

Instruments			
2.904.20	S2.904.20	Drill bit for quick coupling Ø 1.5 mm; spiral 30 mm; 2-flute; length 105 mm	1
2.977.06	–	Drill guide for angle-stable screwing 2.0; for drills Ø 1.5 mm; length 45 mm	2
2.9401.15A	–	Hex screwdriver with handle for screws Ø 1.5 mm and Ø 2.0 mm; wrench width 1.5 mm; detachable	1
2.953.51	–	Gauge with clasp for screws Ø 2.0 mm; measuring range 44 mm	1
2.954.01	–	Screw forceps; self holding	1
2.948.15A	S2.948.15A	Screwdriver for extraction for quick coupling; wrench size 1.5 mm	1
10.428.01	–	Bending forceps for mini plates 1.5 mm – 2.0 mm	1

Wires			
6.031.08	S6.031.08	Kirschner wire with trocar point and round end; Ø 0.8 mm x 150 mm; stainless steel	5


Container			
19.507.00	–	Perforated autoclavable container with inset for instr./impl. angle stable plate MF	1

Optional			
10.305.05	–	T-handle with quick coupling; length 93.5 mm	–
2.9355.35/20	–	Silicone handle with ratchet and AO coupling; length 165 mm; Ø 20 mm	–
2.9345.15	–	Torque screwdriver handle; 0.7 Nm	–

Set MF plate as IH T SET / IMPLANTS sterile

angle-stable; for screws \varnothing 2.0 mm; titanium

non-sterile packaging Set no. 19.508.00 / sterile packaging Set no. 19.508.00PS+



 Internal hexagon
Set no. 19.508.00
non-sterile packaging