

Olecranon Plate 3.5 spike

standard; with spikes; plate thickness 2.0 mm; for screws Ø 3.5 mm; titanium












Item no. 5.874.04



Item no. 5.874.06



Item no. 5.874.10

Item no.	Item description	EAN
5.874.04	olecranon plate with spikes 3.5; 4 holes; 68 x 10 mm; spikes radius 6; titanium	 4047075010837
5.874.05	olecranon plate with spikes 3.5; 5 holes; 80 x 10 mm; spikes radius 6; titanium	 4047075010844
5.874.06	olecranon plate with spikes 3.5; 6 holes; 92 x 10 mm; spikes radius 6; titanium	 4047075010912
5.874.07	olecranon plate with spikes 3.5; 7 holes; 104 x 10 mm; spikes radius 6; titanium	 4047075010929
5.874.08	olecranon plate with spikes 3.5; 8 holes; 116 x 10 mm; spikes radius 6; titanium	 4047075010943
5.874.09	olecranon plate with spikes 3.5; 9 holes; 128 x 10 mm; spikes radius 6; titanium	 4047075010714
5.874.10	olecranon plate with spikes 3.5; 10 holes; 140 x 10 mm; spikes radius 6; titanium	 4047075010721
5.874.12	olecranon plate with spikes 3.5; 12 holes; 164 x 10 mm; spikes radius 6; titanium	 4047075010981
5.874.14	olecranon plate with spikes 3.5; 14 holes; 188 x 10 mm; spikes radius 6; titanium	 4047075010998