


Set Patella Plate 3.5 vas

Variable angle stable; Titanium; Set no. 19.671.010



Item no. 5.8513.11





Item no. 19.671.010  Internal hexagon

Item no.	Item description	Quant.
Implants - plates		
5.8513.11	Patella plate with conical thread \varnothing 3.5 mm; head \varnothing 5.0 mm; 11 holes; preshaped; variable angle stable; titanium	2

Instruments		
2.9405.25	hex screwdriver for \varnothing 2.7 mm - \varnothing 4.0 mm; wrench size 2.5 mm; conical; length 215 mm; silicon	1
2.977.12	drill guide with handle for variable angle-stable screwing \varnothing 3.5 mm; with head \varnothing 5.0 mm	1
2.977.13	drill guide for angle stable screws \varnothing 3.5 mm with head \varnothing 5.0 mm; length 44 mm	2
2.904.06	drill bit for quick coupling; \varnothing 2.5 mm x 125 mm; spiral length 60 mm; double spiral	1
2.953.60	gauge with clasp for screws \varnothing 3.5 mm; \varnothing 4.0 mm; measuring range 60 mm	1
2.954.01	screw forceps; self holding	1


Wire		
6.031.10	Kirschner wire with trocar point and round end; \varnothing 1.0 mm x 150 mm; stainless steel	5

Container		
19.670.010	perforated autoclavable container with inset for instruments and implants patella plate; variable angle stable	1

Item no.	Item description	Quant.
Implants - screws		
Cortical screws with conical head thread \varnothing 3.5 mm; fully threaded; head \varnothing 5.0 mm; self-tapping; titanium  Internal hexagon 		
3.136.12	length 12 mm	5
3.136.14	length 14 mm	5
3.136.16	length 16 mm	5
3.136.18	length 18 mm	5
3.136.20	length 20 mm	5
3.136.22	length 22 mm	5
3.136.24	length 24 mm	5
3.136.26	length 26 mm	5
3.136.28	length 28 mm	5
3.136.30	length 30 mm	5
3.136.32	length 32 mm	5
3.136.36	length 36 mm	5
3.136.40	length 40 mm	5

Set Patella Plate 3.5 vas

Variable angle stable; Titanium; Set no. 19.671.010



Item no. 19.671.010

