

Distal radius plate 2.4 vas Watershed

variable angle-stable; for screws Ø 2.4 mm; titanium



Item no. 5.08373.03







Item no. 5.08373.04



Item no. 5.08473.03



Item no. 5.08473.04

Item no.	Item description	EAN
5.08373.03	dist.radius watershed 2.4; 7/3 holes; 23x50mm; narrow; vas; left	 4 047075 189489
5.08373.04	dist.radius watershed 2.4; 7/4 holes; 23x59mm; narrow; vas; left	 4 047075 189496
5.08473.03	dist.radius watershed 2.4; 7/3 holes; 23x50mm; narrow; vas; right	 4 047075 189526
5.08473.04	dist.radius watershed 2.4; 7/4 holes; 23x59mm; narrow; vas; right	 4 047075 189533

Information

You can use the plate with the instruments of the following sets: 19.117.00 Set Distal Radius Plate 2.4 vas and 19.024.001 Set Instruments 2.4



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498-560
Fax +49 36738 498-569

Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com

AnSp | Edition 01 | 2020-07
All pictured products are CE compliant.