

# Distal Radius Plate 2.4 vas

variable angle-stable; for screws  $\varnothing$  2.4 mm; Titanium



Item no. 5.08173.06
























Item no. 5.08193.06



Item no. 5.08273.06



Item no. 5.08293.06

Item no.	Item description	EAN
5.08173.02	dist.radius plate 2.4; 7/2 holes; 23x46mm; narrow; vas; left	 4 047075 178063
5.08173.03	dist.radius plate 2.4; 7/3 holes; 23x55mm; narrow; vas; left	 4 047075 178070
5.08173.04	dist.radius plate 2.4; 7/4 holes; 23x63mm; narrow; vas; left	 4 047075 178087
5.08173.06	dist.radius plate 2.4; 7/6 holes; 23x80mm; narrow; vas; left	 4 047075 178094
5.08173.08	dist.radius plate 2.4; 7/8 holes; 23x97mm; narrow; vas; left	 4 047075 178100
5.08173.10	dist.radius plate 2.4; 7/10 holes; 23x114mm; narrow; vas; left	 4 047075 178117
5.08173.12	dist.radius plate 2.4; 7/12 holes; 23x131mm; narrow; vas; left	 4 047075 178124
5.08273.02	dist.radius plate 2.4; 7/2 holes; 23x46mm; narrow; vas; right	 4 047075 178209
5.08273.03	dist.radius plate 2.4; 7/3 holes; 23x55mm; narrow; vas; right	 4 047075 178216
5.08273.04	dist.radius plate 2.4; 7/4 holes; 23x63mm; narrow; vas; right	 4 047075 178223
5.08273.06	dist.radius plate 2.4; 7/6 holes; 23x80mm; narrow; vas; right	 4 047075 178230
5.08273.08	dist.radius plate 2.4; 7/8 holes; 23x97mm; narrow; vas; right	 4 047075 178247
5.08273.10	dist.radius plate 2.4; 7/10 holes; 23x114mm; narrow; vas; right	 4 047075 178254
5.08273.12	dist.radius plate 2.4; 7/12 holes; 23x131mm; narrow; vas; right	 4 047075 178261
5.08193.02	dist.radius plate 2.4; 9/2 holes; 28x46mm; broad; vas; left	 4 047075 178131
5.08193.03	dist.radius plate 2.4; 9/3 holes; 28x55mm; broad; vas; left	 4 047075 178148
5.08193.04	dist.radius plate 2.4; 9/4 holes; 28x63mm; broad; vas; left	 4 047075 178155
5.08193.06	dist.radius plate 2.4; 9/6 holes; 28x80mm; broad; vas; left	 4 047075 178162
5.08193.08	dist.radius plate 2.4; 9/8 holes; 28x97mm; broad; vas; left	 4 047075 178179
5.08193.10	dist.radius plate 2.4; 9/10 holes; 28x114mm; broad; vas; left	 4 047075 178186
5.08193.12	dist.radius plate 2.4; 9/12 holes; 28x131mm; broad; vas; left	 4 047075 178193

# Distal Radius Plate 2.4 vas

variable angle-stable; for screws  $\varnothing$  2.4 mm; Titanium

Item no.	Item description	EAN
5.08293.02	dist.radius plate 2.4; 9/2 holes; 28x46mm; broad; vas; right	 4 047075 178278
5.08293.03	dist.radius plate 2.4; 9/3 holes; 28x55mm; broad; vas; right	 4 047075 178285
5.08293.04	dist.radius plate 2.4; 9/4 holes; 28x63mm; broad; vas; right	 4 047075 178292
5.08293.06	dist.radius plate 2.4; 9/6 holes; 28x80mm; broad; vas; right	 4 047075 178308
5.08293.08	dist.radius plate 2.4; 9/8 holes; 28x97mm; broad; vas; right	 4 047075 178315
5.08293.10	dist.radius plate 2.4; 9/10 holes; 28x114mm; broad; vas; right	 4 047075 178322
5.08293.12	dist.radius plate 2.4; 9/12 holes; 28x131mm; broad; vas; right	 4 047075 178339