

# Variable Angle-stable Distal Radius Plates 2.7

for Screws Ø 2.7 mm; Titanium



Item no. 5.9823.54



Item no. 5.9833.54



Item no. 5.9843.74



Item no. 5.9853.74

Item no.	Item description
----------	------------------

## Variable angle-stable distal radius plates - head width 23 mm

5.9823.52	Distal radius plate with conical thread Ø 2.7 mm; 5/ 2 hole; 43 x 10 mm; right; narrow; vas; titanium
5.9823.53	Distal radius plate with conical thread Ø 2.7 mm; 5/ 3 hole; 51 x 10 mm; right; narrow; vas; titanium
5.9823.54	Distal radius plate with conical thread Ø 2.7 mm; 5/ 4 hole; 60 x 10 mm; right; narrow; vas; titanium
5.9823.56	Distal radius plate with conical thread Ø 2.7 mm; 5/ 6 hole; 77 x 10 mm; right; narrow; vas; titanium
5.9823.58	Distal radius plate with conical thread Ø 2.7 mm; 5/ 8 hole; 94 x 10 mm; right; narrow; vas; titanium
5.9823.510	Distal radius plate with conical thread Ø 2.7 mm; 5/10 hole; 111 x 10 mm; right; narrow; vas; titanium
5.9823.513	Distal radius plate with conical thread Ø 2.7 mm; 5/13 hole; 136 x 10 mm; right; narrow; vas; titanium
5.9833.52	Distal radius plate with conical thread Ø 2.7 mm; 5/ 2 hole; 43 x 10 mm; left; narrow; vas; titanium
5.9833.53	Distal radius plate with conical thread Ø 2.7 mm; 5/ 3 hole; 51 x 10 mm; left; narrow; vas; titanium
5.9833.54	Distal radius plate with conical thread Ø 2.7 mm; 5/ 4 hole; 60 x 10 mm; left; narrow; vas; titanium
5.9833.56	Distal radius plate with conical thread Ø 2.7 mm; 5/ 6 hole; 77 x 10 mm; left; narrow; vas; titanium
5.9833.58	Distal radius plate with conical thread Ø 2.7 mm; 5/ 8 hole; 94 x 10 mm; left; narrow; vas; titanium
5.9833.510	Distal radius plate with conical thread Ø 2.7 mm; 5/10 hole; 111 x 10 mm; left; narrow; vas; titanium
5.9833.513	Distal radius plate with conical thread Ø 2.7 mm; 5/13 hole; 136 x 10 mm; left; narrow; vas; titanium

## Variable angle-stable distal radius plates - head width 27 mm

5.9843.72	Distal radius plate with conical thread Ø 2.7 mm; 7/ 2 hole; 43 x 10 mm; right; broad; vas; titanium
5.9843.73	Distal radius plate with conical thread Ø 2.7 mm; 7/ 3 hole; 51 x 10 mm; right; broad; vas; titanium
5.9843.74	Distal radius plate with conical thread Ø 2.7 mm; 7/ 4 hole; 60 x 10 mm; right; broad; vas; titanium
5.9843.76	Distal radius plate with conical thread Ø 2.7 mm; 7/ 6 hole; 77 x 10 mm; right; broad; vas; titanium
5.9843.78	Distal radius plate with conical thread Ø 2.7 mm; 7/ 8 hole; 94 x 10 mm; right; broad; vas; titanium
5.9843.710	Distal radius plate with conical thread Ø 2.7 mm; 7/10 hole; 111 x 10 mm; right; broad; vas; titanium
5.9843.713	Distal radius plate with conical thread Ø 2.7 mm; 7/13 hole; 136 x 10 mm; right; broad; vas; titanium
5.9853.72	Distal radius plate with conical thread Ø 2.7 mm; 7/ 2 hole; 43 x 10 mm; left; broad; vas; titanium
5.9853.73	Distal radius plate with conical thread Ø 2.7 mm; 7/ 3 hole; 51 x 10 mm; left; broad; vas; titanium
5.9853.74	Distal radius plate with conical thread Ø 2.7 mm; 7/ 4 hole; 60 x 10 mm; left; broad; vas; titanium
5.9853.76	Distal radius plate with conical thread Ø 2.7 mm; 7/ 6 hole; 77 x 10 mm; left; broad; vas; titanium
5.9853.78	Distal radius plate with conical thread Ø 2.7 mm; 7/ 8 hole; 94 x 10 mm; left; broad; vas; titanium
5.9853.710	Distal radius plate with conical thread Ø 2.7 mm; 7/10 hole; 111 x 10 mm; left; broad; vas; titanium
5.9853.713	Distal radius plate with conical thread Ø 2.7 mm; 7/13 hole; 136 x 10 mm; left; broad; vas; titanium