

# Distal radius plate 3.5 Twin







angle-stable; for screws  $\varnothing$  3.5 mm; titanium



Item no. 5.876.53

Item no.	Item description	EAN
----------	------------------	-----

## Distal radius plate 3.5 Twin as

5.876.53	Dist. radius plate "Twin" 3.5; 5/3 holes; 51x10mm; R; thickness 3,0mm head/1,7mm shaft	 4047075010882
5.876.55	Dist. radius plate "Twin" 3.5; 5/5 holes; 68x10mm; R; thickness 3,0mm head/1,7mm shaft	 4047075011032
5.876.513	Dist. radius plate "Twin" 3.5; 5/13 holes; 136x10mm; R; thickness 3,0mm head/1,7mm shaft	 4047075090945
5.877.53	Dist. radius plate "Twin" 3.5; 5/3 holes; 51x10mm; L; thickness 3,0mm head/1,7mm shaft	 4047075011216
5.877.54	Dist. radius plate "Twin" 3.5; 5/4 holes; 59x10mm; L; thickness 3,0mm head/1,7mm shaft	 4047075011223
5.877.55	Dist. radius plate "Twin" 3.5; 5/5 holes; 68x10mm; L; thickness 3,0mm head/1,7mm shaft	 4047075011230