

# Distal radius plate 3.5

angle-stable; for screws Ø 3.5 mm; titanium



Item no. 5.876.34



Item no. 5.882.34

Item no.	Item description	EAN
----------	------------------	-----

### Distal radius plate 3.5 Classic as

5.876.33	Dist. radius plate Classic 3.5; 3/3 holes; 51x10mm; R; thickness 3,0mm head/1,7mm shaft	4047075011018
5.876.34	Dist. radius plate Classic 3.5; 3/4 holes; 59x10mm; R; thickness 3,0mm head/1,7mm shaft	4047075011025
5.876.35	Dist. radius plate Classic 3.5; 3/5 holes; 68x10mm; R; thickness 3,0mm head/1,7mm shaft	4047075010783
5.876.36	Dist. radius plate Classic 3.5; 3/6 holes; 76x10mm; R; thickness 3,0mm head/1,7mm shaft	4047075011056
5.876.38	Dist. radius plate Classic 3.5; 3/8 holes; 93x10mm; R; thickness 3,0mm head/1,7mm shaft	4047075011070
5.876.310	Dist. radius plate Classic 3.5; 3/10 holes; 110x10mm; R; thickness 3,0mm head/1,7mm shaft	4047075001989
5.876.313	Dist. radius plate Classic 3.5; 3/13 holes; 136x10mm; R; thickness 3,0mm head/1,7mm shaft	4047075001965
5.877.33	Dist. radius plate Classic 3.5; 3/3 holes; 51x10mm; L; thickness 3,0mm head/1,7mm shaft	4047075011162
5.877.34	Dist. radius plate Classic 3.5; 3/4 holes; 59x10mm; L; thickness 3,0mm head/1,7mm shaft	4047075011179
5.877.36	Dist. radius plate Classic 3.5; 3/6 holes; 76x10mm; L; thickness 3,0mm head/1,7mm shaft	4047075010905
5.877.38	Dist. radius plate Classic 3.5; 3/8 holes; 93x10mm; L; thickness 3,0mm head/1,7mm shaft	4047075011193
5.877.310	Dist. radius plate Classic 3.5; 3/10 holes; 110x10mm; L; thickness 3,0mm head/1,7mm shaft	4047075001996
5.877.313	Dist. radius plate Classic 3.5; 3/13 holes; 136x10mm; L; thickness 3,0mm head/1,7mm shaft	4047075001972

### Distal radius plate 3.5 Minor as

5.882.33	Dist. radius plate Minor 3.5; 3/3 holes; 51x10mm; R; thickness 3,0mm head/1,5mm shaft	4047075061600
5.882.34	Dist. radius plate Minor 3.5; 3/4 holes; 59x10mm; R; thickness 3,0mm head/1,5mm shaft	4047075061617
5.882.36/2	Dist. radius plate Minor 3.5/3.5; 3/6 holes; 76x10mm; R; thickness 3,0mm head/2,0mm shaft	4047075111916
5.882.38/2	Dist. radius plate Minor 3.5/3.5; 3/8 holes; 93x10mm; R; thickness 3,0mm head/2,0mm shaft	4047075111923
5.883.33	Dist. radius plate Minor 3.5; 3/3 holes; 51x10mm; L; thickness 3,0mm head/1,5mm shaft	4047075061624
5.883.34	Dist. radius plate Minor 3.5; 3/4 holes; 59x10mm; L; thickness 3,0mm head/1,5mm shaft	4047075061631
5.883.36/2	Dist. radius plate Minor 3.5/3.5; 3/6 holes; 76x10mm; L; thickness 3,0mm head/2,0mm shaft	4047075111930
5.883.38/3	Dist. radius plate Minor 3.5/3.5; 3/8 holes; 93x10mm; L; thickness 3,0mm head/3,0mm shaft	4047075122233

# Distal radius plate 3.5

angle-stable; for screws Ø 3.5 mm; titanium



Item no. 5.876.44















Item no.	Item description	EAN
<b>Distal radius plate 3.5 Quadro as</b>		
5.876.43	Dist. radius plate Quadro 3.5; 4/3 holes; 51x10mm; R; thickness 3,0mm head/1,7mm shaft	 4047075032198
5.876.43/2	Dist. radius plate Quadro 3.5; 4/3 holes; 51x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075057412
5.876.44	Dist. radius plate Quadro 3.5; 4/4 holes; 59x10mm; R; thickness 3,0mm head/1,7mm shaft	 4047075032204
5.876.44/2	Dist. radius plate Quadro 3.5; 4/4 holes; 59x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075057429
5.876.45/2	Dist. radius plate Quadro 3.5; 4/5 holes; 68x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075057887
5.876.46/2	Dist. radius plate Quadro 3.5; 4/6 holes; 76x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075111879
5.876.48/2	Dist. radius plate Quadro 3.5; 4/8 holes; 93x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075111886
5.876.48/3	Dist. radius plate Quadro 3.5; 4/8 holes; 93x10mm; R; thickness 3,0mm head/3,0mm shaft	 4047075056040
5.877.43	Dist. radius plate Quadro 3.5; 4/3 holes; 51x10mm; L; thickness 3,0mm head/1,7mm shaft	 4047075032167
5.877.44	Dist. radius plate Quadro 3.5; 4/4 holes; 59x10mm; L; thickness 3,0mm head/1,7mm shaft	 4047075032174
5.877.44/2	Dist. radius plate Quadro 3.5; 4/4 holes; 59x10mm; L; thickness 3,0mm head/2,0mm shaft	 4047075057498
5.877.45/2	Dist. radius plate Quadro 3.5; 4/5 holes; 68x10mm; L; thickness 3,0mm head/2,0mm shaft	 4047075057870
5.877.46/2	Dist. radius plate Quadro 3.5; 4/6 holes; 76x10mm; L; thickness 3,0mm head/2,0mm shaft	 4047075111893
5.877.48/2	Dist. radius plate Quadro 3.5; 4/8 holes; 93x10mm; L; thickness 3,0mm head/2,0mm shaft	 4047075111909
5.877.48/3	Dist. radius plate Quadro 3.5; 4/8 holes; 93x10mm; L; thickness 3,0mm head/3,0mm shaft	 4047075056057

# Distal radius plate 3.5

angle-stable; for screws Ø 3.5 mm; titanium



Item no. 5.886.34

Item no.	Item description	EAN
<b>Distal radius plate 3.5 Vario as</b>		
5.886.310/3	Dist. radius plate Vario 3.5/3.5; 3/10 holes; 110x10mm; R; thickness 3,0mm head/3,0mm shaft	 4047075111206
5.886.313/3	Dist. radius plate Vario 3.5/3.5; 3/13 holes; 136x10mm; R; thickness 3,0mm head/3,0mm shaft	 4047075111213
5.886.33	Dist. radius plate Vario 3.5/3.5; 3/3 holes; 51x10mm; R; thickness 3,0mm head/1,7mm shaft	 4047075011254
5.886.34	Dist. radius plate Vario 3.5/3.5; 3/4 holes; 59x10mm; R; thickness 3,0mm head/1,7mm shaft	 4047075011315
5.886.36/2	Dist. radius plate Vario 3.5/3.5; 3/6 holes; 76x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075111183
5.886.38/2	Dist. radius plate Vario 3.5/3.5; 3/8 holes; 93x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075111190
5.886.38/3	Dist. radius plate Vario 3.5/3.5; 3/8 holes; 93x10mm; R; thickness 3,0mm head/3,0mm shaft	 4047075122240
5.887.310/3	Dist. radius plate Vario 3.5/3.5; 3/10 holes; 110x10mm; L; thickness 3,0mm head/3,0mm shaft	 4047075111244
5.887.313/3	Dist. radius plate Vario 3.5/3.5; 3/13 holes; 136x10mm; L; thickness 3,0mm head/3,0mm shaft	 4047075111251
5.887.33	Dist. radius plate Vario 3.5/3.5; 3/3 holes; 51x10mm; L; thickness 3,0mm head/1,7mm shaft	 4047075011322
5.887.34	Dist. radius plate Vario 3.5/3.5; 3/4 holes; 59x10mm; L; thickness 3,0mm head/1,7mm shaft	 4047075011339
5.887.36/2	Dist. radius plate Vario 3.5/3.5; 3/6 holes; 76x10mm; L; thickness 3,0mm head/2,0mm shaft	 4047075111220
5.887.38/2	Dist. radius plate Vario 3.5/3.5; 3/8 holes; 93x10mm; L; thickness 3,0mm head/2,0mm shaft	 4047075111237
5.887.38/3	Dist. radius plate Vario 3.5/3.5; 3/8 holes; 93x10mm; L; thickness 3,0mm head/1,7mm shaft	 4047075122257