

Distal radius plate 3.5 vas
























variable angle-stable; for screws Ø 3.5 mm; titanium



Item no. 5.8763.34



Item no. 5.8823.34

Item no.	Item description	EAN
Distal radius plate 3.5 Classic vas		
5.8763.33	Distal radius plate Classic 3.5; 3/3 holes; 51x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075111428
5.8763.34	Distal radius plate Classic 3.5; 3/4 holes; 59x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075111435
5.8763.36	Distal radius plate Classic 3.5; 3/6 holes; 76x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075121106
5.8763.38	Distal radius plate Classic 3.5; 3/8 holes; 93x10mm; R; thickness 3,0mm head/3,0mm shaft	 4047075121113
5.8763.310	Distal radius plate Classic 3.5; 3/10 holes; 110x10mm; R; thickness 3,0mm head/3,0mm shaft	 4047075121083
5.8763.313	Distal radius plate Classic 3.5; 3/13 holes; 136x10mm; R; thickness 3,0mm head/3,0mm shaft	 4047075121090
5.8773.33	Distal radius plate Classic 3.5; 3/3 holes; 51x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075111442
5.8773.34	Distal radius plate Classic 3.5; 3/4 holes; 59x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075111459
5.8773.36	Distal radius plate Classic 3.5; 3/6 holes; 76x10mm; L; thickness 3,0mm head/2,0mm shaft	 4047075121144
5.8773.38	Distal radius plate Classic 3.5; 3/8 holes; 93x10mm; L; thickness 3,0mm head/3,0mm shaft	 4047075121151
5.8773.310	Distal radius plate Classic 3.5; 3/10 holes; 110x10mm; L; thickness 3,0mm head/3,0mm shaft	 4047075121120
5.8773.313	Distal radius plate Classic 3.5; 3/13 holes; 136x10mm; L; thickness 3,0mm head/3,0mm shaft	 4047075121137
Distal radius plate 3.5 Minor vas		
5.8823.33	Distal radius plate Minor 3.5; 3/3 holes; 51x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075116492
5.8823.34	Distal radius plate Minor 3.5; 3/4 holes; 59x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075116508
5.8823.36	Distal radius plate Minor 3.5; 3/6 holes; 76x10mm; R; thickness 3,0mm head/2,0mm shaft	 4047075120338
5.8823.38	Distal radius plate Minor 3.5; 3/8 holes; 93x10mm; R; thickness 3,0mm head/3,0mm shaft	 4047075120345
5.8823.310	Distal radius plate Minor 3.5; 3/10 holes; 110x10mm; R; thickness 3,0mm head/3,0mm shaft	 4047075120314
5.8823.313	Distal radius plate Minor 3.5; 3/13 holes; 136x10mm; R; thickness 3,0mm head/3,0mm shaft	 4047075120321
5.8833.33	Distal radius plate Minor 3.5; 3/3 holes; 51x10mm; L; thickness 3,0mm head/2,0mm shaft	 4047075116515
5.8833.34	Distal radius plate Minor 3.5; 3/4 holes; 59x10mm; L; thickness 3,0mm head/2,0mm shaft	 4047075116522
5.8833.36	Distal radius plate Minor 3.5; 3/6 holes; 76x10mm; L; thickness 3,0mm head/2,0mm shaft	 4047075120376
5.8833.38	Distal radius plate Minor 3.5; 3/8 holes; 93x10mm; L; thickness 3,0mm head/3,0mm shaft	 4047075120383
5.8833.310	Distal radius plate Minor 3.5; 3/10 holes; 110x10mm; L; thickness 3,0mm head/3,0mm shaft	 4047075120352
5.8833.313	Distal radius plate Minor 3.5; 3/13 holes; 136x10mm; L; thickness 3,0mm head/3,0mm shaft	 4047075120369