

Set Distal Radius plates 3.5 vas


for screws Ø 3.5 mm; titanium; Set no. 19.1123.00



Item no. 5.8763.33 Classic



Item no. 5.8823.33 Minor



Set no. 19.1123.00  Internal hexagon

Choose 4 pieces in total between ...

Item no.	Item description	Quant.
Implants – plates		
Distal Radius plate 3.5 Classic vas, titanium		
5.8763.33	3/3 holes; 51x10mm; vas; R; titanium; thickness 3,0mm head/2,0mm shaft	1
5.8763.34	3/4 holes; 59x10mm; vas; R; titanium; thickness 3,0mm head/2,0mm shaft	1
5.8773.33	3/3 holes; 51x10mm; vas; L; titanium; thickness 3,0mm head/2,0mm shaft	1
5.8773.34	3/4 holes; 59x10mm; vas; L; titanium; thickness 3,0mm head/2,0mm shaft	1

and/or:

Implants – plates		
Distal Radius plate 3.5 Minor vas, titanium		
5.8823.33	3/3 holes; 51 x 10mm; vas; R; titanium; thickness 3,0mm head/2,0mm shaft	1
5.8823.34	3/4 holes; 59 x 10mm; vas; R; titanium; thickness 3,0mm head/2,0mm shaft	1
5.8833.33	3/3 holes; 51 x 10mm; vas; L; titanium; thickness 3,0mm head/2,0mm shaft	1
5.8833.34	3/4 holes; 59 x 10mm; vas; L; titanium; thickness 3,0mm head/2,0mm shaft	1

Implants – screws		
Cortical screw Ø 3.5 mm; fully threaded; self tapping; titanium		
3.132.14	length 14 mm	5
3.132.16	length 16 mm	5
3.132.18	length 18 mm	5
3.132.20	length 20 mm	5
3.132.22	length 22 mm	5
3.132.24	length 24 mm	5

Implants – screws		
Angle-stable Cortical screw Ø 3.5 mm; fully threaded; self tapping; titanium		
3.133.16	length 16 mm	3
3.133.18	length 18 mm	3
3.133.20	length 20 mm	3
3.133.22	length 22 mm	3
3.133.24	length 24 mm	3
3.133.26	length 26 mm	3
3.133.28	length 28 mm	3
3.133.30	length 30 mm	3

Set Distal Radius plates 3.5 vas


for screws Ø 3.5 mm; titanium; Set no. 19.1123.00



Item no. 3.132.XX



Item no. 3.133.XX



Set no. 19.1123.00  Internal hexagon


Item no.	Item description	Quant.
Instruments		
10.500.27A	drill guide with handle for small fragment, handle length 55mm	1
2.904.05	drill bit for quick coupling, Ø2.5x105 mm, spiral length 50mm, double spiral	1
2.9406.25	hex screwdriver with handle for screws Ø2.7mm - Ø4.0 mm, AF 2.5 mm, conical	1
2.953.60	gauge with clasp for screws Ø3.5 mm, Ø4.0 mm, measuring range 60mm	1
2.954.01	screw forceps, self holding	1
2.977.01	drill guide for angle stable screwing, length 45mm, small fragment	1
2.977.08	drill guide w.handle for variable angle stable screwing, small fragment, length 45 mm	1

Item no.	Item description	Quant.
Wires		
6.031.18	Kirschner wire with trocar point and round end, Ø1.8x150 mm, st.st.	5

Container		
19.113.00	perforated autoclavable container with inset for instruments/implants radius-set	1

Set Distal Radius plates 3.5 vas

for screws Ø 3.5 mm; titanium; Set no. 19.1123.00



Set no. 19.1123.00