

Reco plates 3.5

standard; for screws Ø 3.5 mm; stainless steel



Item no. 4.312.06



Item no. 4.313.06



Item no. 4.328.06

Item no.	Item description	EAN
----------	------------------	-----

Reco plate 3.5 – plate thickness 2.6 mm

4.312.03	reco-plate 3.5; 3 holes; 34 x 10 mm; straight; stainless steel	
4.312.04	reco-plate 3.5; 4 holes; 46 x 10 mm; straight; stainless steel	
4.312.05	reco-plate 3.5; 5 holes; 58 x 10 mm; straight; stainless steel	
4.312.06	reco-plate 3.5; 6 holes; 70 x 10 mm; straight; stainless steel	
4.312.07	reco-plate 3.5; 7 holes; 82 x 10 mm; straight; stainless steel	
4.312.08	reco-plate 3.5; 8 holes; 94 x 10 mm; straight; stainless steel	
4.312.09	reco-plate 3.5; 9 holes; 106 x 10 mm; straight; stainless steel	
4.312.10	reco-plate 3.5; 10 holes; 118 x 10 mm; straight; stainless steel	
4.312.11	reco-plate 3.5; 11 holes; 130 x 10 mm; straight; stainless steel	
4.312.12	reco-plate 3.5; 12 holes; 142 x 10 mm; straight; stainless steel	
4.312.13	reco-plate 3.5; 13 holes; 154 x 10 mm; straight; stainless steel	
4.312.14	reco-plate 3.5; 14 holes; 166 x 10 mm; straight; stainless steel	
4.312.15	reco-plate 3.5; 15 holes; 178 x 10 mm; straight; stainless steel	
4.312.16	reco-plate 3.5; 16 holes; 190 x 10 mm; straight; stainless steel	
4.312.17	reco-plate 3.5; 17 holes; 202 x 10 mm; straight; stainless steel	
4.312.18	reco-plate 3.5; 18 holes; 214 x 10 mm; straight; stainless steel	
4.312.19	reco-plate 3.5; 19 holes; 226 x 10 mm; straight; stainless steel	
4.312.20	reco-plate 3.5; 20 holes; 238 x 10 mm; straight; stainless steel	
4.312.22	reco-plate 3.5; 22 holes; 262 x 10 mm; straight; stainless steel	

Reco plate 3.5 strong – plate thickness 3.2 mm

4.313.04	reco-plate 3.5; 4 holes; 46 x 10 mm; straight; strong; stainless steel	
4.313.05	reco-plate 3.5; 5 holes; 58 x 10 mm; straight; strong; stainless steel	
4.313.06	reco-plate 3.5; 6 holes; 70 x 10 mm; straight; strong; stainless steel	
4.313.07	reco-plate 3.5; 7 holes; 82 x 10 mm; straight; strong; stainless steel	
4.313.08	reco-plate 3.5; 8 holes; 94 x 10 mm; straight; strong; stainless steel	



Königsee Implantate GmbH
 OT Aschau | Am Sand 4
 07426 Allendorf | Germany

Fon +49 36738 498-560
 Fax +49 36738 498-569

Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com

KrMo | Edition 02 | 2020-09
 All pictured products are CE compliant.

Reco plates 3.5

standard; for screws Ø 3.5 mm; stainless steel

Item no.	Item description	EAN
4.313.09	reco-plate 3.5; 9 holes; 106 x 10 mm; straight; strong; stainless steel	 4 047075 103867
4.313.10	reco-plate 3.5; 10 holes; 118 x 10 mm; straight; strong; stainless steel	 4 047075 103874
4.313.11	reco-plate 3.5; 11 holes; 130 x 10 mm; straight; strong; stainless steel	 4 047075 118311
4.313.12	reco-plate 3.5; 12 holes; 142 x 10 mm; straight; strong; stainless steel	 4 047075 104048
4.313.13	reco-plate 3.5; 13 holes; 154 x 10 mm; straight; strong; stainless steel	 4 047075 104055
4.313.14	reco-plate 3.5; 14 holes; 166 x 10 mm; straight; strong; stainless steel	 4 047075 103966
4.313.15	reco-plate 3.5; 15 holes; 178 x 10 mm; straight; strong; stainless steel	 4 047075 118328
4.313.16	reco-plate 3.5; 16 holes; 190 x 10 mm; straight; strong; stainless steel	 4 047075 104079
4.313.17	reco-plate 3.5; 17 holes; 202 x 10 mm; straight; strong; stainless steel	 4 047075 118335
4.313.18	reco-plate 3.5; 18 holes; 214 x 10 mm; straight; strong; stainless steel	 4 047075 103973
4.313.19	reco-plate 3.5; 19 holes; 226 x 10 mm; straight; strong; stainless steel	 4 047075 118342
4.313.20	reco-plate 3.5; 20 holes; 238 x 10 mm; straight; strong; stainless steel	 4 047075 103980
4.313.22	reco-plate 3.5; 22 holes; 262 x 10 mm; straight; strong; stainless steel	 4 047075 103997

Reco plate 3.5 bent – plate thickness 2.6 mm

4.328.03	reco-plate 3.5; 3 holes; breadth 10; arcuate; stainless steel	 4 047075 104192
4.328.04	reco-plate 3.5; 4 holes; breadth 10; arcuate; stainless steel	 4 047075 104369
4.328.05	reco-plate 3.5; 5 holes; breadth 10; arcuate; stainless steel	 4 047075 104376
4.328.06	reco-plate 3.5; 6 holes; breadth 10; arcuate; stainless steel	 4 047075 104277
4.328.07	reco-plate 3.5; 7 holes; breadth 10; arcuate; stainless steel	 4 047075 104383
4.328.08	reco-plate 3.5; 8 holes; breadth 10; arcuate; stainless steel	 4 047075 104284
4.328.09	reco-plate 3.5; 9 holes; breadth 10; arcuate; stainless steel	 4 047075 104291
4.328.10	reco-plate 3.5; 10 holes; breadth 10; arcuate; stainless steel	 4 047075 104307
4.328.12	reco-plate 3.5; 12 holes; breadth 10; arcuate; stainless steel	 4 047075 104222
4.328.14	reco-plate 3.5; 14 holes; breadth 10; arcuate; stainless steel	 4 047075 104499
4.328.16	reco-plate 3.5; 16 holes; breadth 10; arcuate; stainless steel	 4 047075 104390
4.328.18	reco-plate 3.5; 18 holes; breadth 10; arcuate; stainless steel	 4 047075 104406

Distal humerus plate Y 3.5 Reco – plate thickness 2.2 mm

4.861.03	Y-reco plate for distal humerus 3.5; 3 holes; breadth 10 mm; stainless steel	 4 047075 048632
4.861.06	Y-reco plate for distal humerus 3.5; 6 holes; breadth 10 mm; stainless steel	 4 047075 048717
4.861.07	Y-reco plate for distal humerus 3.5; 7 holes; breadth 10 mm; stainless steel	 4 047075 037872
4.861.08	Y-reco plate for distal humerus 3.5; 8 holes; breadth 10 mm; stainless steel	 4 047075 048793
4.861.10	Y-reco plate for distal humerus 3.5; 10 holes; breadth 10 mm; stainless steel	 4 047075 048625