

Reco plate 4.5 acetabulum DC

standard; for screws Ø 4.5 mm; stainless steel



Item no. 4.338.06

Item no.	Item description	EAN
Reco plate 4.5 acetabulum DC – plate thickness 3.2 mm		
4.338.03	reco-plate for acetabulum with DC holes 4.5; 3 holes; 46 x 11.5 mm; stainless steel	4 047075 104635
4.338.04	reco-plate for acetabulum with DC holes 4.5; 4 holes; 62 x 11.5 mm; stainless steel	4 047075 104642
4.338.05	reco-plate for acetabulum with DC holes 4.5; 5 holes; 78 x 11.5 mm; stainless steel	4 047075 104666
4.338.06	reco-plate for acetabulum with DC holes 4.5; 6 holes; 94 x 11.5 mm; stainless steel	4 047075 104673
4.338.07	reco-plate for acetabulum with DC holes 4.5; 7 holes; 110 x 11.5 mm; stainless steel	4 047075 104871
4.338.08	reco-plate for acetabulum with DC holes 4.5; 8 holes; 126 x 11.5 mm; stainless steel	4 047075 104888
4.338.09	reco-plate for acetabulum with DC holes 4.5; 9 holes; 142 x 11.5 mm; stainless steel	4 047075 104789
4.338.10	reco-plate for acetabulum with DC holes 4.5; 10 holes; 158 x 11.5 mm; stainless steel	4 047075 104901
4.338.11	reco-plate for acetabulum with DC holes 4.5; 11 holes; 174 x 11.5 mm; stainless steel	4 047075 104796
4.338.12	reco-plate for acetabulum with DC holes 4.5; 12 holes; 190 x 11.5 mm; stainless steel	4 047075 104802
4.338.13	reco-plate for acetabulum with DC holes 4.5; 13 holes; 206 x 11.5 mm; stainless steel	4 047075 104819
4.338.14	reco-plate for acetabulum with DC holes 4.5; 14 holes; 222 x 11.5 mm; stainless steel	4 047075 104703
4.338.15	reco-plate for acetabulum with DC holes 4.5; 15 holes; 238 x 11.5 mm; stainless steel	4 047075 104994
4.338.16	reco-plate for acetabulum with DC holes 4.5; 16 holes; 254 x 11.5 mm; stainless steel	4 047075 104918