

# Set Angle-stable Lateral T-Tibia Head Plate

for Screws Ø 4.5/6.0 mm; Large Fragment; Titanium; Set no. 19.408.00



Item no. 5.777.44



Item no. 5.778.44



Set no. 19.408.00

Internal hexagon


Item no.	Item description	Quant.
<b>Implants - plates</b>		
<b>T-tibia head plate with conical thread 6.0/4.5; titanium; left</b>		
5.777.44	4 holes; 95 mm x 14 mm	1
5.777.46	6 holes; 127 mm x 14 mm	1
5.777.48	8 holes; 159 mm x 14 mm	1
5.777.410	10 holes; 191 mm x 14 mm	1


<b>T-tibia head plate with conical thread 6.0/4.5; titanium; right</b>		
5.778.44	4 holes; 95 mm x 14 mm	1
5.778.46	6 holes; 127 mm x 14 mm	1
5.778.48	8 holes; 159 mm x 14 mm	1
5.778.410	10 holes; 191 mm x 14 mm	1

<b>Instruments</b>		
2.9406.35	hex screwdriver with handle for screws Ø 4.5 mm - Ø 7.0 mm; AF 3.5 mm; length 215 mm, PPSU	1
2.977.02	drill guide for angle-stable screwing; length 70 mm; large fragment	2
2.904.09	drill bit for quick coupling; Ø 3.2 mm x 195 mm; spiral length 60 mm; double spiral	1
2.954.01	screw forceps; self holding	1
2.953.90	gauge with clasp for screws Ø 4.5 mm; Ø 6.0 mm; measuring range 90 mm	1

<b>Wire</b>		
6.031.18	Kirschner wire with trocar point and round end; Ø 1.8 mm x 150 mm; stainless steel	5


<b>Container</b>		
19.400.00/8	perforated autoclavable container with inset for instruments and proximal tibia head plates	1

Item no.	Item description	Quant.
<b>Implants - screws</b>		
<b>Cortical screw with conical head thread Ø 4.5 mm; fully threaded self tapping; titanium</b>		
		
Internal hexagon		
3.545.18	length 18 mm	5
3.545.20	length 20 mm	5
3.545.22	length 22 mm	5
3.545.24	length 24 mm	5
3.545.26	length 26 mm	5
3.545.28	length 28 mm	5
3.545.30	length 30 mm	5
3.545.32	length 32 mm	5
3.545.34	length 34 mm	5
3.545.36	length 36 mm	5
3.545.38	length 38 mm	5
3.545.40	length 40 mm	5
3.545.42	length 42 mm	5
3.545.45	length 45 mm	5
3.545.50	length 50 mm	5
3.545.55	length 55 mm	5



<b>Cancelous screw with conical head thread Ø 6.0 mm; fully threaded self tapping; titanium</b>		
		
Internal hexagon		
3.163.55	length 55 mm	4
3.163.60	length 60 mm	4
3.163.65	length 65 mm	4
3.163.70	length 70 mm	4
3.163.75	length 75 mm	4
3.163.80	length 80 mm	4
3.163.85	length 85 mm	4
3.163.90	length 90 mm	4

# Set Angle-stable Lateral T-Tibia Head Plate

for Screws  $\varnothing$  4.5/6.0 mm; Large Fragment; Titanium; Set no. 19.408.00



Set no. 19.408.00  Internal hexagon

Item no.	Item description	Quant.
<b>Implants - Screws</b>		
Cortical screw $\varnothing$ 4.5 mm fully threaded; self tapping titanium		  Internal hexagon
3.152.18	length 18 mm	5
3.152.20	length 20 mm	5
3.152.22	length 22 mm	5
3.152.24	length 24 mm	5
3.152.26	length 26 mm	5
3.152.28	length 28 mm	5
3.152.30	length 30 mm	5
3.152.32	length 32 mm	5
3.152.34	length 34 mm	5
3.152.36	length 36 mm	5
3.152.38	length 38 mm	5
3.152.40	length 40 mm	5
3.152.42	length 42 mm	5
3.152.45	length 45 mm	5
3.152.50	length 50 mm	5