

Set TIS™ Screws


with instruments; TiOB; titanium; Set no. 19.681.00



Item no. 3.531.90



Item no. 3.530.90



Set no. 19.681.00  Internal hexagon

Item no.	Item description	Quant.
Implants - screws		
TIS™ Screws; full thread; Ø 7.5 mm, inner Ø 2.9 mm; self-drilling; self-tapping; titanium		
3.530.50	Length 50 mm	2
3.530.55	Length 55 mm	2
3.530.60	Length 60 mm	2
3.530.65	Length 65 mm	2
3.530.70	Length 70 mm	2
3.530.75	Length 75 mm	2
3.530.80	Length 80 mm	2
3.530.85	Length 85 mm	2
3.530.90	Length 90 mm	2
3.530.95	Length 95 mm	2
3.530.100	Length 100 mm	2
3.530.105	Length 105 mm	2
3.530.110	Length 110 mm	2
3.530.115	Length 115 mm	2
3.530.120	Length 120 mm	2
3.530.125	Length 125 mm	2
3.530.130	Length 130 mm	2

Instruments		
2.908.50	drill bit w.coupling, Ø5x300, spiral length 55, cannulated Ø2.9, four-spiral	1
2.9435.75	hex screwdriver with handle for Ø7.5, cannulated, IØ2.9, wrench size 4, L 324mm	1
2.952.200	gauge for cannulated screws Ø7.5, measuring range 200, length 275mm	1
2.954.01	screw forceps, self holding	1

Item no.	Item description	Quant.
Implants - screws		
TIS™ Screws; thread 32mm; Ø 7.5 mm, inner Ø 2.9 mm; self-drilling; self-tapping; back-tapping; titanium		
3.531.50	Length 50 mm	2
3.531.55	Length 55 mm	2
3.531.60	Length 60 mm	2
3.531.65	Length 65 mm	2
3.531.70	Length 70 mm	2
3.531.75	Length 75 mm	2
3.531.80	Length 80 mm	2
3.531.85	Length 85 mm	2
3.531.90	Length 90 mm	2
3.531.95	Length 95 mm	2
3.531.100	Length 100 mm	2
3.531.105	Length 105 mm	2
3.531.110	Length 110 mm	2
3.531.115	Length 115 mm	2
3.531.120	Length 120 mm	2
3.531.125	Length 125 mm	2
3.531.130	Length 130 mm	2

Implants – washer		
3.200.04	washer, outside Ø16 / inside Ø6.7, ti	5

Wires		
2.980.28	guide wire, Ø2.8x300, thread 5, st. st.	5



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498-560
Fax +49 36738 498-569

Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com

Set TIS™ Screws


with instruments; TiOB; titanium; Set no. 19.681.00

Container

19.680.00

perforated autoclavable container
w.inset f.instr.a.TIS-screws Ø7,5

1



Set no. 19.681.00

Information:

The name TIS™ derives from the German designation "Transiliosakralschraube", which translates into English as "transiliosacral screw". Products with the add-on "™" are trademarks of Königsee Implantate.

TiOB® is a surface treatment that brings optimized properties.



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498-560
Fax +49 36738 498-569

Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com

AnSp | Edition 04 | 2021-04
All pictured products are CE compliant.