

# Set Angle-stable Medial Tibia Head Plate


for Screws Ø 4.5/6.0 mm; Titanium; Set no. 19.409.00



Item no. 5.7192.36



Item no. 5.7182.36



Set no. 19.409.00



 Internal hexagon



Item no.	Item description	Quant.
<b>Implants - plates</b>		
<b>Medial tibia head plate Ø 4.5 mm - 6.0 mm; titanium; angle-stable</b>		
5.7192.34	4 holes; 77 mm x 14 mm; left	1
5.7192.36	6 holes; 109 mm x 14 mm; left	1
5.7192.38	8 holes; 141 mm x 14 mm; left	1
5.7182.34	4 holes; 77 mm x 14 mm; right	1
5.7182.36	6 holes; 109 mm x 14 mm; right	1
5.7182.38	8 holes; 141 mm x 14 mm; right	1

<b>Instruments</b>		
2.9405.35	hex screwdriver for Ø 4.5 mm - Ø 7.0 mm; wrench size 3.5 mm; length 215 mm; silicone	1
2.904.09	drill bit for quick coupling; Ø 3.2 mm x 195 mm; spiral length 60 mm; double spiral	1
2.977.02	drill guide for angle stable screwing; length 70 mm; big fragment	2
2.953.90	gauge with clasp for screws Ø 4.5 mm; Ø 6.0 mm; measuring range 90 mm	1
2.954.01	screw forceps; self holding	1
10.502.32	universal drill guide; Ø 3.2 mm	1

<b>Wire</b>		
6.031.18	Kirschner wire with trocar point and round end; Ø 1.8 mm x 150 mm; stainless steel	5


<b>Container</b>		
19.400.10	perforated autoclavable container with inset for instruments and medial tibia head plates	1

Item no.	Item description	Quant.
<b>Implants - screws</b>		
<b>Canellous screw with conical head thread Ø 6.0 mm; fully threaded; self tapping; titanium</b>		
 Internal hexagon		
3.163.55	length 55 mm	4
3.163.60	length 60 mm	4
3.163.65	length 65 mm	4
3.163.70	length 70 mm	4
3.163.75	length 75 mm	4

<b>Cortical screw with conical head thread Ø 4.5 mm; fully threaded; self tapping; titanium</b>		
		
 Internal hexagon		
3.545.18	length 18 mm	5
3.545.20	length 20 mm	5
3.545.22	length 22 mm	5
3.545.24	length 24 mm	5
3.545.26	length 26 mm	5
3.545.28	length 28 mm	5
3.545.30	length 30 mm	5
3.545.32	length 32 mm	5
3.545.34	length 34 mm	5
3.545.36	length 36 mm	5
3.545.38	length 38 mm	5
3.545.40	length 40 mm	5
3.545.42	length 42 mm	5
3.545.45	length 45 mm	5
3.545.50	length 50 mm	5
3.545.55	length 55 mm	5


# Set Angle-stable Medial Tibia Head Plate

for Screws Ø 4.5/6.0 mm; Titanium; Set no. 19.409.00




Set no. 19.409.00



 Internal hexagon



Cortical screws



Cancellous screws

Item no.	Item description	Quant.
Implants - screws		
Cortical screw Ø 4.5 mm; fully threaded; self tapping; titanium		  Internal hexagon
3.152.18	length 18 mm	3
3.152.20	length 20 mm	3
3.152.22	length 22 mm	3
3.152.24	length 24 mm	3
3.152.26	length 26 mm	3
3.152.28	length 28 mm	3
3.152.30	length 30 mm	3
3.152.32	length 32 mm	3
3.152.34	length 34 mm	3
3.152.36	length 36 mm	3
3.152.38	length 38 mm	3
3.152.40	length 40 mm	3
3.152.42	length 42 mm	3
3.152.45	length 45 mm	3
3.152.50	length 50 mm	3