

Set Distal ulna plate 2.0


angle-stable; titanium; Set no. 19.506.00



Item no. 5.888.43







Item no. 5.889.43



Set no. 19.506.00  Internal hexagon

Item no.	Item description	Quant.
Implants – plates*		
distal ulna plate with conical thread 2.0; titanium		
5.888.43	4/3 holes; 38 x 6 mm; right	1
5.888.46	4/6 holes; 60 x 6 mm; right	1
5.889.43	4/3 holes; 38 x 6 mm; left	1
5.889.46	4/6 holes; 60 x 6 mm; left	1

Implants – screws		
cortical screw with conical head thread; Ø 2.0 mm; fully threaded; self tapping; titanium		
		
 Internal hexagon		
3.114.10	length 10 mm	5
3.114.12	length 12 mm	5
3.114.14	length 14 mm	5
3.114.16	length 16 mm	5

cortical screw; Ø 2.0 mm; fully threaded; self tapping; titanium		
		
 Internal hexagon		
3.110.10	length 10 mm	5
3.110.12	length 12 mm	5
3.110.14	length 14 mm	5


Item no.	Item description	Quant.
Instruments		
2.904.20	drill bit for quick coupling; Ø 1.5 x 105; spiral length 30; double spiral	1
2.977.06	drill guide for angle stable screwing 2.0; length 45	2
2.9401.15A	hex screwdriver with handle for screw Ø 1.5 and Ø 2.0; wrench width 1.5; detachable	1
2.953.51	gauge with clasp for screws Ø 2.0; measuring range 44	1
2.954.01	screw forceps; self holding	1

Container		
19.507.00	perforated autoclavable container w.inset f.instr./impl. angle stable plate MF	1

*Optionally available: 5.888.45 und 5.889.45

Set Distal ulna plate 2.0

angle-stable; titanium; Set no. 19.506.00



Set no. 9.506.00  Internal hexagon