

Distal ulna plate 2.0

angle-stable; for screws \varnothing 2.0 mm; plate thickness 1.5 mm; titanium



Item no. 5.888.43



Item no. 5.889.43

Item no.	Item description	EAN
5.888.43	distal ulna plate with conical thread 2.0; 4/3 holes; 38 x 6 mm; right; titanium	 4047075116669
5.888.45	distal ulna plate with conical thread 2.0; 4/5 holes; 53 x 6 mm; right; titanium	 4047075110735
5.888.46	distal ulna plate with conical thread 2.0; 4/6 holes; 60 x 6 mm; right; titanium	 4047075116676
5.889.43	distal ulna plate with conical thread 2.0; 4/3 holes; 38 x 6 mm; left; titanium	 4047075116645
5.889.45	distal ulna plate with conical thread 2.0; 4/5 holes; 53 x 6 mm; left; titanium	 4047075110728
5.889.46	distal ulna plate with conical thread 2.0; 4/6 holes; 60 x 6 mm; left; titanium	 4047075116652