

VaWiKo® EPI plate


variable angle-stable; for screws Ø 4.0 mm; titanium; TiOB®; plate thickness 1.9 mm

non-sterile packaging / sterile packaging



with spikes

Item no. 5.6983.2122
Item no. S5.6983.2122



without spikes

Item no. 5.6983.2160
Item no. S5.6983.2160

Item no.	EAN	Item description	Item no. sterile	EAN sterile
VaWiKo® EPI plate				
5.6983.2122		VaWiKo® EPI plate for hemiepiphysodesis; hole spacing 12; 2 spikes	S5.6983.2122	
5.6983.2162		VaWiKo® EPI plate for hemiepiphysodesis; hole spacing 16; 2 spikes	S5.6983.2162	
5.6983.2120		VaWiKo® EPI plate for hemiepiphysodesis; hole spacing 12; without spikes	S5.6983.2120	
5.6983.2160		VaWiKo® EPI plate for hemiepiphysodesis; hole spacing 16; without spikes	S5.6983.2160	

Information

The system name VaWiKo® is derived from the German designation „Variabel winkelstabile Kombination“, which translates into English as „variable angle-stable combination“. This system combines conventional plating technologies with angle-stable locking in variable angles.

VaWiKo® is a registered trademark of Königsee Implantate GmbH.

TiOB® is a surface treatment that brings optimized properties.



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498-560
Fax +49 36738 498-569

export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com

KrMo | Edition 04 | 2026-01
All pictured products are CE compliant.