


Set VaWiKo® Lateral distal fibula plate 3.5

variable angle-stable; titanium; Set no. 19.474.00/2







Item no. 5.07813.09



Set no. 19.474.00/2  Internal hexagon

Item no.	Item description	Quant.
Implants – plates		
VaWiKo® lateral distal fibula plate 3.5; titanium		
5.07813.06	4/2 holes; length 83 mm	1
5.07813.07	4/3 holes; length 96 mm	1
5.07813.09	4/5 holes; length 122 mm	1

Implants – screws		
Cortical screw with conical head thread Ø 3.5 mm; fully threaded; self tapping; head 5.0; titanium		  Internal hexagon
3.136.10	length 10 mm	4
3.136.12	length 12 mm	4
3.136.14	length 14 mm	4
3.136.16	length 16 mm	4
3.136.18	length 18 mm	4
3.136.20	length 20 mm	4
3.136.22	length 22 mm	4
3.136.24	length 24 mm	4

Item no.	Item description	Quant.
Implants – screws		
Cortical screw Ø 3.5 mm self tapping, fully threaded; titanium		  Internal hexagon
3.138.10	length 10 mm	2
3.138.12	length 12 mm	2
3.138.14	length 14 mm	2
3.138.16	length 16 mm	2
3.138.18	length 18 mm	2
3.138.20	length 20 mm	2
3.138.22	length 22 mm	2
3.138.24	length 24 mm	2
3.138.36	length 36 mm	2
3.138.38	length 38 mm	2
3.138.40	length 40 mm	2
3.138.42	length 42 mm	2
3.138.44	length 44 mm	2
3.138.46	length 46 mm	2
3.138.48	length 48 mm	2
3.138.50	length 50 mm	2

Set VaWiKo® Lateral distal fibula plate 3.5


variable angle-stable; titanium; Set no. 19.474.00/2


Item no.	Item description	Quant.
Instruments		
2.904.06	drill bit for quick coupling, Ø 2.5 x 125 mm, spiral length 60 mm, double spiral	1
2.904.10	drill bit for quick coupling, Ø 3.5 x 110 mm, spiral length 40 mm, double spiral	1
2.977.12	drill guide with handle for variable angle stable screwing 3.5 with head Ø 5.0 mm	1
2.977.13	drill guide for angle stable screws 3.5 with head Ø 5.0 mm, length 44 mm	2
10.520.30	double drill guide Ø 2.5 mm/ Ø 3.5 mm, spring bushing	1
2.9406.25	hex screwdriver with handle for screws Ø 2.7 mm - Ø 4.0 mm, AF 2.5, conical	1
2.953.60	gauge with clasp for screws Ø 3.5 mm, Ø 4.0 mm, measuring range 60 mm	1
2.924.35	tap for quick coupling for cortical screws, Ø 3.5 x 110 mm	1
2.954.01	screw forceps, self holding	1

Item no.	Item description	Quant.
Accessories		
3.200.02	washer, outside Ø 7.0 mm / inside Ø 3.8 mm, titanium	3

Wire		
6.031.16	Kirschner wire with trocar point and round end, Ø 1.6 x 150 mm, stainless steel	5

Container		
19.475.00/2	perforated autoclav. container with inset for instruments/ implants vas, plate and screws up to 50 mm in length	1



Set no. 19.474.00/2  Internal hexagon