

EN VaWiKo[®] proximal Humerus plate 4.0/3.5

variable angle-stable; for screws Ø 4.0/3.5 mm; titanium



Item no. 5.07913.84

Item no.	Item description	EAN
VaWiKo[®] proximal Humerus plate 4.0/3.5; titanium		
5.07913.83	VaWiKo [®] proximal Humerus plate 4.0/3.5; 8/3 hole; length 77 mm	4 047075 167029
5.07913.84	VaWiKo [®] proximal Humerus plate 4.0/3.5; 8/4 hole; length 93 mm	4 047075 167036
5.07913.85	VaWiKo [®] proximal Humerus plate 4.0/3.5; 8/5 hole; length 109 mm	4 047075 167043
5.07913.86	VaWiKo [®] proximal Humerus plate 4.0/3.5; 8/6 hole; length 125 mm	4 047075 168231
5.07913.87	VaWiKo [®] proximal Humerus plate 4.0/3.5; 8/6 hole; length 141 mm	4 047075 167050
5.07913.88	VaWiKo [®] proximal Humerus plate 4.0/3.5; 8/8 hole; length 157 mm	4 047075 167067

Information

The system name VaWiKo[®] is derived from the German designation „Variabel winkelstabile Kombination“, which translates into English as „variable angle-stable combination“. This system combines conventional plating technologies with angle-stable locking in variable angles.

VaWiKo[®] is a registered trademark of Königsee Implantate GmbH.



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498-560
Fax +49 36738 498-569

Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com

HeBo | Edition 02 | 2019-10
All pictured products are CE compliant.