

WiKo® Anterior Clavicula Plate 3.5

medial and lateral; angle-stable; for screws Ø 3.5 mm; titanium; TiOB®



Item no. 5.3802.06



Item no. 5.3812.08

Item no.	Item description
----------	------------------

WiKo® anterior medial clavicula plate 3.5; titanium; TiOB®

5.3802.06	WiKo® anterior medial clavicula plate 3.5; 6 hole; length 80 mm
5.3802.07	WiKo® anterior medial clavicula plate 3.5; 7 hole; length 91 mm
5.3802.08	WiKo® anterior medial clavicula plate 3.5; 8 hole; length 103 mm

WiKo® anterior lateral clavicula plate 3.5; titanium; TiOB®

5.3812.08	WiKo® anterior lateral clavicula plate 3.5; 8 hole; length 105 mm
5.3812.09	WiKo® anterior lateral clavicula plate 3.5; 9 hole; length 117 mm
5.3812.10	WiKo® anterior lateral clavicula plate 3.5; 10 hole; length 128 mm

Information:

The system name WiKo® is derived from the German designation "Winkelstabile Kombination", which translates into English as "angle stable combination". This is a plating technology which offers stable locking at a specified angle.

Products which feature the "®" symbol are registered trademarks of Königsee Implantate GmbH.