

WiKo® Proximal Femur-osteotomy plate 3.5/6.0

DVO Pediatric Hip System; children plate; stainless steel; Set no. 19.561.10



Item no. 4.6402.0404




Item no. 4.6402.0803




Item no. 4.6402.0603






Item no. 4.6402.0403



Set no. 19.561.10  Internal hexagon


Item no.	Item description	Quant.
Implants – plates		
WiKo® proximal femur-osteotomy plate; stainless steel		
4.6402.0403	3.5/6.0; 3 holes; 70 mm; arc depth 4	1
4.6402.0404	3.5/6.0; 4 holes; 81 mm; arc depth 4	1
4.6402.0603	3.5/6.0; 3 holes; 73 mm; arc depth 6	1
4.6402.0604	3.5/6.0; 4 holes; 84 mm; arc depth 6	1
4.6402.0803	3.5/6.0; 3 holes; 77 mm; arc depth 8	1

Item no.	Item description	Quant.
Implants – screws		
Cortical screw Ø 3.5 mm; fully threaded; self tapping; angle-stable; stainless steel		
2.133.14	length 14 mm	3
2.133.16	length 16 mm	3
2.133.18	length 18 mm	3
2.133.20	length 20 mm	3
2.133.22	length 22 mm	3
2.133.24	length 24 mm	3
2.133.26	length 26 mm	3
2.133.28	length 28 mm	3
2.133.30	length 30 mm	3
2.133.32	length 32 mm	3
2.133.34	length 34 mm	3
2.133.36	length 36 mm	3
2.133.38	length 38 mm	3

Item no.	Item description	Quant.
Implants – screws		
Cannulated screw Ø 6.0 mm; fully thread.; self drilling; self tapping; inside-Ø 2.2 mm; angle-stable; stainl. steel		
2.540.28	length 28 mm	2
2.540.30	length 30 mm	2
2.540.34	length 34 mm	2
2.540.38	length 38 mm	2
2.540.42	length 42 mm	2
2.540.46	length 46 mm	2
2.540.50	length 50 mm	2
2.540.52	length 52 mm	2
2.540.56	length 56 mm	2
2.540.60	length 60 mm	2
2.540.65	length 65 mm	2
Cortical screw Ø 3.5 mm; fully threaded; self tapping; stainless steel		
2.132.14	length 14 mm	2
2.132.16	length 16 mm	2
2.132.18	length 18 mm	2
2.132.20	length 20 mm	2
2.132.22	length 22 mm	2
2.132.24	length 24 mm	2
2.132.26	length 26 mm	2
2.132.28	length 28 mm	2
2.132.30	length 30 mm	2
2.132.32	length 32 mm	2
2.132.34	length 34 mm	2
2.132.36	length 36 mm	2
2.132.38	length 38 mm	2

WiKo® Proximal Femur-osteotomy plate 3.5/6.0


DVO Pediatric Hip System; children plate; stainless steel; Set no. 19.561.10



upper level, femur-osteotomy plates and triangular positioning plate



upper level, Instruments Large fragment



lower level, Instruments Small fragment



lower level, targeting device and guide sleeve

Item no.	Item description	Quant.
Instruments		
10.230.01	triangular positioning plate 80°/70°/30°	1
10.230.02	triangular positioning plate 90°/50°/40°	1
10.230.03	triangular positioning plate 100°/60°/20°	1
10.230.20	targeting device; variable	1
10.230.21	holding instrument for targeting device	1
10.230.24	guide sleeve for screwdriver; Ø 15 mm x 80 mm; inside Ø 7.0 mm	1
10.230.22	cleaning adapter for targeting device 10.230.20; Ø 6.2 mm with hose connection	1
2.954.01	screw forceps; self holding	1
Large fragment		
2.942.02	hex screwdriver with T-handle for screws Ø 4.5 - 7.0 mm; cannulated; inside Ø 2.2 mm	1
2.908.40	drill bit with triangular shank; Ø 4.0 x 230 mm; spiral length 40 mm; cannulated; inside Ø 2.1 mm; triple spiral	1
2.952.230	gauge for cannulated screws Ø 6.5 mm and Ø 7.0 mm; length 230 mm	1
Small fragment		
2.9405.25	hex screwdriver for screws Ø 2.7 - 4.0 mm; wrench width 2.5; conical; length 215 mm	1
2.904.22	drill bit for quick coupling; Ø 2.5 x 155 mm; spiral length 60 mm; double spiral	1
2.977.04	drill guide for angle stable screwing; length 68 mm	2
10.500.04	drill guide for neutral and compression position; Ø 2.5 mm	1
2.953.70	gauge with clasp for screws with conical head thread Ø 3.5 mm and Ø 4.0 mm; measuring range 70 mm	1

Item no.	Item description	Quant.
Wire		
6.031.20	Kirschner wire with trocar point and round end; Ø 2.0 mm x 150 mm; stainless steel	5
6.041.20/230	Kirschner wire with trocar point and round end; Ø 2.0 mm x 230 mm; stainless steel	5
Container		
19.560.00	perforated autoclavable container with inset for instruments for WiKo®-DVO; Pediatric Hip System	1

Information:


The system name WiKo® derives from the description „angle-stable combination“. This is a plating technology that provides a stable lock at a given angle.

WiKo® is a registered trademark of Königsee Implantate GmbH in Germany.


DVO - Derotation-Varisation-Osteotomy

WiKo® Proximal Femur-osteotomy plate 3.5/6.0

DVO Pediatric Hip System; children plate; stainless steel; Set no. 19.561.10



Set no. 19.561.10, upper level  Internal hexagon



Set no. 19.561.10, lower level  Internal hexagon