

Set Hand Arthrodesis Plate vas

variable angle-stable; for screws Ø 2.7/3.5 mm; head 5.0; titanium; Set no. 19.671.020



Item no. 5.9963.3340




Item no. 5.9963.3342



Item no. 5.9963.3350





Item no. 5.9963.4540





Set no. 19.671.020 Internal hexagon

Item no.	Item description	Quant.
Implants – plates		
Hand Arthrodesis Plate vas; variable angle-stable; titanium		
5.9963.3340	3/3 holes; 4°; offset 0	1
5.9963.3342	3/3 holes; 4°; offset 2	1
5.9963.3350	3/3 holes; 5°; offset 0	1
5.9963.4540	4/5 holes; 4°; offset 0	1

Implants – screws		
Cortical screw; Ø 2.7; fully threaded; self tapping; flat head; titanium		
 Internal hexagon		
3.128.12	length 12 mm	2
3.128.14	length 14 mm	2
3.128.16	length 16 mm	2
3.128.18	length 18 mm	2
3.128.20	length 20 mm	2
3.128.22	length 22 mm	2
3.128.24	length 24 mm	2
3.128.26	length 26 mm	2
3.128.28	length 28 mm	2
3.128.30	length 30 mm	2

Cortical screw with conical head thread; Ø 2.7; fully threaded; s.t.; titanium		
 Internal hexagon		
3.125.12	length 12 mm	5
3.125.14	length 14 mm	5
3.125.16	length 16 mm	5
3.125.18	length 18 mm	5

Item no.	Item description	Quant.
Implants – screws		
cortical screw; Ø 3.5; fully threaded; head 5.0; s.t.; titanium; TiOB®		
 Internal hexagon		
3.138.12	length 12 mm	2
3.138.14	length 14 mm	2
3.138.16	length 16 mm	2
3.138.18	length 18 mm	2
3.138.20	length 20 mm	2
3.138.22	length 22 mm	2
3.138.24	length 24 mm	2
3.138.26	length 26 mm	2
3.138.28	length 28 mm	2
3.138.30	length 30 mm	2
3.138.32	length 32 mm	2
3.138.36	length 36 mm	2
3.138.40	length 40 mm	2

cortical screw with conical head thread; Ø 3.5; fully threaded; head 5.0; s.t.; titanium; TiOB®		
 Internal hexagon		
3.136.12	length 12 mm	5
3.136.14	length 14 mm	5
3.136.16	length 16 mm	5
3.136.18	length 18 mm	5
3.136.20	length 20 mm	5
3.136.22	length 22 mm	5
3.136.24	length 24 mm	5
3.136.26	length 26 mm	5

Set Hand Arthrodesis Plate vas

variable angle-stable; for screws Ø 2.7/3.5 mm; head 5.0; titanium; Set no. 19.671.020

Item no.	Item description	Quant.
3.125.20	length 20 mm	5
3.125.22	length 22 mm	5
3.125.24	length 24 mm	5
3.125.26	length 26 mm	5
3.125.28	length 28 mm	5
3.125.30	length 30 mm	5

Wires		
6.031.10	Kirschner wire with trocar point and round end, Ø1x150, st. st.	5


Container		
19.670.020	perforated autoclavable cont.w.inset f.instr./impl.hand arthrodesis plates,vas	1

Item no.	Item description	Quant.
3.136.28	length 28 mm	5
3.136.30	length 30 mm	5
3.136.32	length 32 mm	5
3.136.36	length 36 mm	5
3.136.40	length 40 mm	5

Instruments		
2.9406.25	hex screwdriver with handle for screws Ø2.7-Ø4.0, AF2.5, conical	1
2.977.10	drill guide w.handle for variable angle stable screwing, mini fragment, length 45	1
2.977.11	drill guide for angle stable screwing 2.7, length 44	2
2.977.12	drill guide with handle for variable angle stable screwing 3.5 with head Ø 5.0	1
2.977.13	drill guide for angle stable screws 3.5 with head Ø 5.0, length 44	2
2.904.04	drill bit for quick coupling, Ø2x125, spiral length 30, double spiral	1
2.904.06	drill bit for quick coupling, Ø2.5x125, spiral length 60, double spiral	1
2.953.65	gauge w. clasp for screws Ø2.7, measuring range 60	1
2.954.01	screw forceps, self holding	1

Set Hand Arthrodesis Plate vas

variable angle-stable; for screws \varnothing 2.7/3.5 mm; head 5.0; titanium; Set no. 19.671.020



Set no. 19.671.020  Internal hexagon