

L-AD plate 3.5 skew

Variable angle-stable; plate thickness 2.3 mm; for screws
Ø 3.5 mm; head 5.0 mm; skew; titanium




2/2 holes
33 x 8 mm
left/right



2/3 holes
43 x 8 mm
left/right

LCF plate 3.5

Variable angle-stable; plate thickness 2.3 mm; for screws
Ø 3.5 mm; head 5.0 mm; 8 holes; length 74 mm;
titanium



left/right



Contact details

If you are interested in our
HINDFOOT PLATES,
please feel free to contact us



Königsee Implantate GmbH

OT Aschau | Am Sand 4
07426 Allendorf
Germany

export@koenigsee-implantate.de
+ 49 36738 498-560


or your local sales representatives.

www.koenigsee-implantate.de

Copyright © Königsee Implantate GmbH/DaTi | Edition 03 | 2020-07

HINDFOOT PLATES 3.5

titanium plates range



...especially for 3.5 and 2.7 screws!

 **Königsee**
Implantate

Foot plate AD 3.5

Variable angle-stable; plate thickness 2.3 mm;
for screws Ø 3.5 mm; head 5.0 mm; titanium



4 holes
length 23 and 28 mm
2 variants



Foot plate AD interposition 3.5


Variable angle-stable; plate thickness 2.3 mm;
for screws Ø 3.5 mm; head 5.0 mm; width 20 mm;
titanium




5 holes
length 28 and 33 mm
2 variants

Multi-AD plate 3.5

Variable angle-stable; plate thickness 2.3 mm;
for screws Ø 3.5 mm; head 5.0 mm; width 20 mm;
titanium



6 holes
length 43 mm
long bridge



6 - 12 holes
length 38 - 83 mm
4 variants

T-AD plate 3.5

Variable angle-stable; plate thickness 2.3 mm;
for screws Ø 3.5 mm; head 5.0 mm; titanium



2/2 holes
33 x 8 mm



2/2 holes
38 x 8 mm

T-AD plate 3.5 DC

Variable angle-stable; 1 DC-holes; plate thickness
2.3 mm; for screws Ø 3.5 mm; head 5.0 mm;
titanium



2/2 holes
33 x 8 mm



2/2 holes
38 x 8 mm



2/3 holes
43 x 8 mm

Information The system name VaWiKo® is derived from the German designation „Variabel winkelstabile Kombination“, which translates into English as „variable angle-stable combination“. This system combines conventional plating technologies with angle-stable locking in variable angles. VaWiKo® is a registered trademark of Königsee Implantate GmbH.

Cortical screws

Cortical screw 3.5 head 5; fully threaded; self-tapping; titanium; TiOB®

Internal hexagon




Cortical screw 3.5 head 5; fully threaded; angle-stable; self-tapping; titanium; TiOB®

Internal hexagon



Information

Set Hindfoot plates 3.5 holeØ5 vas can be combined with Set Forefoot plates 2.7 vas.



Set Hindfoot plates 3.5 vas

VaWiKo® Lapidus Plate 3.5

Variable angle-stable; plate thickness 2.0 mm;
for screws Ø 3.5 mm; head 5.0 mm;
41 x 13 mm; titanium



plantar offset 2
left/right

