

Hand plates Mini 2.0 as

angle-stable; for screws Ø 2.0 mm; plate thickness 1.2 mm; titanium



Item no. 5.9012.05



Item no. 5.9022.04



Item no. 5.9052.04S



Item no. 5.9332.36



Item no. 5.9012.13

















Item no. 5.9032.06



Item no. 5.9312.26



Item no. 5.9342.34


Item no.	Item description	EAN
Straight plate Mini 2.0 as		
5.9012.05	straight plate with conical thread 2.0; 5 holes; 29 x 5 mm; titanium	 4047075129799
5.9012.13	straight plate with conical thread 2.0; 13 holes; 77 x 5 mm; titanium	 4047075129805
Straight Plate Mini 2.0 as round bridge		
5.9022.04	mini bridge plate with conical thread 2.0; 4 holes; 29 x 5 mm; titanium	 4047075148806
K-Plate Mini 2.0 as		
5.9032.06	K-plate with conical thread 2.0; 6 holes; length 37 mm; titanium	 4047075148813
Double bridge plate 2.0 as		
5.9052.04S	double bridge plate with conical thread 2.0; 4 holes; 15 x 11 mm; oblique; titanium	 4047075129812
5.9052.07S	double bridge plate with conical thread 2.0; 7 holes; 25 x 10 mm; oblique; titanium	 4047075148820
5.9052.10	double bridge plate with conical thread 2.0; 10 holes; 45 x 11 mm; straight; titanium	 4047075129829
T-Plate Mini 2.0 as narrow		
5.9312.24	mini T-plate with conical thread 2.0/2.0; 2/4 holes; 30 x 5 mm; titanium	 4047075129836
5.9312.26	mini T-plate with conical thread 2.0/2.0; 2/6 holes; 42 x 5 mm; titanium	 4047075158508
5.9312.28	mini T-plate with conical thread 2.0/2.0; 2/8 holes; 54 x 5 mm; titanium	 4047075158515
T-Plate Mini 2.0 as broad		
5.9332.34	mini T-plate with conical thread 2.0/2.0; 3/4 holes; 30 x 5 mm; straight; titanium	 4047075129850
5.9332.36	mini T-plate with conical thread 2.0/2.0; 3/6 holes; 42 x 5 mm; straight; titanium	 4047075158539
5.9332.38	mini T-plate with conical thread 2.0/2.0; 3/8 holes; 54 x 5 mm; straight; titanium	 4047075158546
T-Plate Mini 2.0 as Δ		
5.9342.34	mini T-plate with conical thread 2.0/2.0; 2/5 holes; 33 x 5 mm; titanium	 4047075148837

Hand plates Mini 2.0 as

angle-stable; for screws \varnothing 2.0 mm; plate thickness 1.2 mm; titanium



Item no. 5.9322.34



Item no. 5.9322.36



Item no. 5.9322.38

Item no.	Item description	EAN
5.9322.34	mini T-plate with conical thread 2.0/2.0; 3/4 holes; 33 x 5 mm; y-form; titanium	 4047075129843
5.9322.36	mini T-plate with conical thread 2.0/2.0; 3/6 holes; 45 x 5 mm; y-form; titanium	 4047075149612
5.9322.38	mini T-plate with conical thread 2.0/2.0; 3/8 holes; 57 x 5 mm; y-form; titanium	 4047075149629