


Set Ring plates vas

variable angle-stable; two levels; titanium; Set no. 19.439.00

LEVEL 1



Item no.	Item description	Quant.
Implants – plates		
Ring plate with conical thread; variable angle-stable; titanium		
5.9953.14	Ring plate 2.0 vas; 7 holes; Ø 14.0 mm	1
5.9953.16	Ring plate 2.7 vas; 5 holes; Ø 16.5 mm	1

Implants – screws		
Cortical screw 2.0; fully threaded; self tapping; titanium; titanium Internal hexagon		
3.110.08	length 8 mm	2
3.110.10	length 10 mm	2
3.110.12	length 12 mm	2
3.110.14	length 14 mm	2
3.110.16	length 16 mm	2
3.110.18	length 18 mm	2
3.110.20	length 20 mm	2

Cortical screw with conical head thread 2.0; fully threaded; self tapping; titanium Internal hexagon		
3.114.08	length 8 mm	5
3.114.10	length 10 mm	5
3.114.12	length 12 mm	5
3.114.14	length 14 mm	5
3.114.16	length 16 mm	5
3.114.18	length 18 mm	5
3.114.20	length 20 mm	5

Item no.	Item description	Quant.
Implants – screws		
Cortical screw with conical head thread 2.7; fully threaded; self tapping; titanium Internal hexagon		
3.125.08	length 8 mm	5
3.125.10	length 10 mm	5
3.125.12	length 12 mm	5
3.125.14	length 14 mm	5
3.125.16	length 16 mm	5
3.125.18	length 18 mm	5
3.125.20	length 20 mm	5

Instruments		
2.954.01	screw forceps, self holding	1
2.977.06	drill guide for angle stable screwing 2.0, length 45	2

LEVEL 2


Instruments		
2.904.03	drill bit for quick coupling, Ø2x105, spiral length 30, double spiral	1
2.904.20	drill bit for quick coupling, Ø1.5x105, spiral length 30, double spiral	1
2.9102.14	cutter 80° for quick coupling, for ring plate Ø 14.0	1
2.9102.16	cutter 100° for quick coupling, for ring plate Ø 16.0	1
2.9406.25	hex screwdriver with handle for screws Ø2.7-Ø4.0, AF2.5, conical	1

Set Ring plates vas



variable angle-stable; two levels; titanium; Set no. 19.439.00

Item no.	Item description	Quant.
Container		
19.438.00	perforated autoclavable container with inset for instruments and ring plates	1

Item no.	Item description	Quant.
2.9401.15A	hex screwdriver with handle for screw \varnothing 1.5 and \varnothing 2.0, wrench width 1.5, detachable	1
2.953.51	gauge with clasp for screws \varnothing 2.0, measuring range 44	1
2.953.65	gauge w. clasp for screws \varnothing 2.7, measuring range 60	1
2.977.10	drill guide w.handle for variable angle stable screwing, mini fragment, length 45	1
2.977.11	drill guide for angle stable screwing 2.7, length 44	2
2.977.14	drill guide for variable angle stable screwing 2.0, length 45	1



Set no. 19.439.00 Complete  Internal hexagon



Set no. 19.439.00 Level 1



Set no. 19.439.00 Level 2