

# Straight plates Mini 2.0 as

angle-stable; for screws  $\varnothing$  2.0 mm; stainless steel











Item no. 4.9002.20



Item no. 4.4555.06



Item no. 4.470.06

Item no.	Item description	EAN
<b>Straight plate Mini 2.0 as – plate thickness 1.5 mm</b>		
4.9002.20	straight mini plate with conical thread 2.0; 20 holes; 119 x 5 mm; stainless steel	
<b>Straight plate Mini 2.0 as bridge – plate thickness 2.0 mm</b>		
4.4555.04	straight plate with conical thread and bridge 2.0; 4 holes; 27 x 5 mm; bridge 12; div.; stainless steel	
4.4555.06	straight plate with conical thread and bridge 2.0; 6 holes; 37 x 5 mm; bridge 12; div.; stainless steel	
4.4555.08	straight plate with conical thread and bridge 2.0; 8 holes; 47 x 5 mm; bridge 12; div.; stainless steel	
4.4555.10	straight plate with conical thread and bridge 2.0; 10 holes; 57 x 5 mm; bridge 12; div.; stainless steel	
<b>Straight plate Mini 2.0 DC as bridge 16 – plate thickness 2.0 mm</b>		
4.470.04	straight plate with conical thread and bridge 2.0; 4 holes; 31 x 5 mm; bridge 16; div.; stainless steel	
4.470.06	straight plate with conical thread and bridge 2.0; 6 holes; 41 x 5 mm; bridge 16; div.; stainless steel	
4.470.08	straight plate with conical thread and bridge 2.0; 8 holes; 51 x 5 mm; bridge 16; div.; stainless steel	
4.470.10	straight plate with conical thread and bridge 2.0; 10 holes; 61 x 5 mm; bridge 16; div.; stainless steel	